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The diagram consists of three main elements: a large light purple oval on the left containing the text '滅菌剤' (Disinfectant); a blue mound in the center containing the text '水分' (Moisture); and a grey rectangular block at the bottom containing the text '器材' (Equipment). A thin purple line connects the top-right corner of the purple oval to the top-left corner of the text box on the right.

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Warning health care facilities that adenosine triphosphate (ATP) test strips should not be used to assess duodenoscope cleaning. To date, the FDA has not evaluated them for effectiveness for assessing duodenoscope reprocessing. Manufacturers of ATP test strips are advised to submit data to support the legal marketing of these strips for this use.

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<https://www.fda.gov/medical-devices/reprocessing-reusable-medical-devices/infections-associated-reprocessed-duodenoscopes>
Infections Associated with Reprocessed Duodenoscopesより引用

In early 2015 the FDA emphasized the importance of training, oversight, and competency testing for reprocessing staff and close attention to optimal cleaning of the elevator mechanism on the leading end of duodenoscopes.

48 In May 2015 an FDA advisory panel emphasized the importance of ensuring availability of duodenoscopes for clinical care and expressed concern that existing

reprocessing practices, even when performed appropriately, appeared to be insufficient for ensuring safety of side-viewing duodenoscopes.⁵⁰ In August 2015 the

FDA went on to advise that together with “strict adherence to the duodenoscope manufacturer’s reprocessing instructions” all centers using duodenoscopes should

“closely evaluate whether they have the expertise, training and resources to implement one or more of several different options,” including microbiologic culturing,

ethylene oxide sterilization, use of a liquid chemical sterilant processing system, and repeated HLD

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Presentation by

William Rutala, PhD, MPH:

University North Carolina Healthcare

Also referenced

Flexible Endoscope Reprocessing and the Importance of AAMI ST91 presentation by Mary Drosnock from Healthmark Industries

Summary by Barb Keller RN BSN CIC

EH Spaulding believed that how an object will be disinfected depended on the object's intended use (developed in 1968)

CRITICAL: objects which enter normally sterile tissue or the vascular system or through which blood flows should be sterile.

SEMICRITICAL: objects that touch mucous membranes or skin that is not intact require a disinfection procedure (high-level disinfection that kills all microorganisms but high numbers of bacterial spores).

NONCRITICAL: objects that touch only intact skin require low-level disinfection (or non-germicidal detergent).

“The Good, The Bad and The Ugly” OUR RESPONSIBILITY TO THE FUTURE

Presentation by William Rutala, PhD, MPH: University North Carolina Healthcare

https://higherlogicdownload.s3.amazonaws.com/APIC/01cbfd6f-d5f9-448b-ac39-3d0525af7514/UploadedImages/Sept%202016%20Meeting/Disinfection_and_Sterilization.pptx

EH Spaulding believed that how an object will be disinfected depended on the object's intended use (now modified)

CRITICAL: objects which **directly or secondarily (i.e., via a mucous membrane such as duodenoscope, cystoscope, bronchoscope)** enter normally sterile tissue or the vascular system or through which blood flows should be sterile.

SEMICRITICAL: objects that touch mucous membranes or skin that is not intact require a disinfection procedure (high-level disinfection that kills all microorganisms but high numbers of bacterial spores).

NONCRITICAL: objects that touch only intact skin require low-level disinfection (or non-germicidal detergent).

“The Good, The Bad and The Ugly” OUR RESPONSIBILITY TO THE FUTURE

Presentation by William Rutala, PhD, MPH: University North Carolina Healthcare

https://higherlogicdownload.s3.amazonaws.com/APIC/01cbfd6f-d5f9-448b-ac39-3d0525af7514/UploadedImages/Sept%202016%20Meeting/Disinfection_and_Sterilization.pptx

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