

Archer® LiquidPlex™ Protocol for Illumina®

Released July 9, 2020

Notices

Limitations of Use

For Research Use Only. Not for use in diagnostic procedures. Not intended to be used for treatment of human or animal diseases.

Safety data sheets pertaining to this product are available upon request.

Safety Notices



Caution symbols denote critical steps in the procedure where risk of assay failure or damage to the product itself could occur if not carefully observed.



Stop symbols indicate where this procedure may be safely suspended and resumed at a later time without risk of compromised assay performance. Make note of these steps and plan your workflow accordingly.



Reminder symbols call attention to minor details that may be easily overlooked and compromise to the procedure resulting in decreased assay performance.

Revision History

Version	Revisions
LA090.2 Released July 9, 2020	Updated formatting and LiquidPlex™ product name.

In This Guide

This protocol is a guide to using Archer NGS library preparation kits, MBC adapters and gene enrichment panels for targeted sequencing of select genes and regions of interest using next-generation sequencing (NGS).

Overview – page 4

This section contains the intended use statement, test principle and a high level overview of the workflow as well as a description of how the required reagents are supplied.

Before Getting Started – page 8

This section contains critical guidance for the successful implementation of the protocol and library preparation. This should be read and understood before laboratory work is initiated.

Protocol – page 10

This section is the step-by-step protocol describing how to perform the workflow.

Additional Resources

View videos and additional resources for Archer products at:

<http://archerdx.com/support>

Technical Support

Visit <http://archerdx.com/support/faqs> for a list of helpful answers to frequently asked questions or contact us directly at tech@archerdx.com

Overview

Intended Use

The Archer LiquidPlex Protocol is intended for use with Archer reagent kits and corresponding target enrichment panels to produce high-complexity libraries for use with Illumina next-generation sequencing (NGS) instruments.

Sequencing data produced by this method should be processed using Archer Analysis software - a integrated bioinformatics suite that leverages Anchored Multiplex PCR (AMP™) chemistry to detect unique sequence fragments, enabling error correction, read deduplication, and high-confidence alignment and variant calling. Archer Analysis begins with demultiplexed FASTQ files from the sequencer and produces both high-level and detailed variant calling. Raw text and BAM file outputs provide full pipeline transparency.

Test Principle

AMP is a patented target enrichment chemistry engineered to produce sequencer-ready NGS libraries from low input, highly fragmented or damaged nucleic acid templates, such as FFPE and cfDNA inputs. Archer lyophilized reagent kits provide a rapid and scalable method for targeted NGS library preparation with single-use reactions in an optimized single-day workflow that supports the use of multichannel pipettors. AMP technology can be used for targeted RNA or DNA sequencing and genotyping applications.


AMP chemistry enriches specific genomic sequences from highly fragmented nucleic acid templates by utilizing a single gene-specific primer (GSP) instead of traditional amplicon-based chemistry where the template must have both opposing GSP binding sites intact for PCR amplification. Unamplified template molecules (cDNA, gDNA or cfDNA) are directly ligated using Archer MBC Adapters to incorporate a sample index, molecule barcode, and a universal priming site to enable single-GSP target enrichment.

Archer Library preparation reagents include:


- Archer ctDNA reagents in lyophilized format for each step of library preparation.
- Gene specific primers (GSPs) that target panel specific regions of interest during PCR amplification.
- Archer Molecular Barcode (MBC) Adapters are proprietary adapters that tag each unique molecule with a barcode and common region prior to amplification.

Modular Assay Format

Archer Kits include library preparation reagents and assay-specific liquid primers, which are used in conjunction with Archer MBC Adapters to construct sequencing-ready libraries from total nucleic acid or DNA samples. See individual product inserts for assay targets and read depth requirements.



Protocol Overview



Supplied Reagents

LiquidPlex Reagents, for Illumina-HS (SK0119)

Store at 2°C to 8°C

Allow pouches to reach room temperature before opening.

Description	Part Number	Quantity
Complete End Repair	SA0205	1 pouch (8 reactions/ 8 tube strip)
Ligation Step 1	SA0196	
Ligation Step 2	SA0197	
First PCR for Illumina-HS	SA0359	
Second PCR Reactions 1 thru 8 for Illumina-HS	SA0360	
Ligation Cleanup Beads	SA0210	1 tube/bottle (sufficient for processing 8 samples)
Ligation Cleanup Buffer	SA0209	
500 mM Tris-HCl, pH 8.0	SA0020	
Ultra-Pure Water	SA0213	
Ultra-Pure Water for Ethanol Dilution	SA0022	

LiquidPlex Frozen Components (part number varies)

Store at -30°C to -10°C

Component Name	Part Number	Quantity
LiquidPlex GSP1	Refer to product insert	Varies (Refer to product insert)
LiquidPlex GSP2	Refer to product insert	

Materials Required But Not Supplied

Additional Materials Required for Archer Library Preparation

Description	Supplier	Part Number
Archer MBC Adapters for Illumina	ArcherDX	Varies
Agencourt® AMPure® XP Beads	Beckman Coulter	A63880
100% Ethanol (ACS grade)	Various	-
RNase AWAY™	Thermo Fisher Scientific	7003
KAPA Universal Library Quantification Kit	KAPA Biosystems	KK4824
MiSeq® or NextSeq® Reagent Kit (300 cycle minimum)	Illumina	-
PhiX Control v3	Illumina	FC-110-3001
200mM Tris pH 7.0 (NextSeq loading)	Various	-
5mM NaOH (ACS grade)	Various	-
Standard PCR Thermal Cycler	Various	-
Real-Time PCR Thermal Cycler	Various	-
qPCR tubes	Various	-
0.2mL PCR tubes	Various	-
DynaMag™-96 Side Magnet	Thermo Fisher Scientific	12331D
Microcentrifuge	Various	-
Plate centrifuge	Various	-
Pipettes (P10, P20, P200 and P1000)	Pipetman P10, P20, P200, P1000 or equivalent	-
Sterile, nuclease-free aerosol barrier pipette tips	Various	-
Vortex Mixer	Various	-
PCR tube cooling block	Various	-
Gloves	Various	-
Qubit® 3.0 Fluorometer	Thermo Fisher Scientific	Q33216
Qubit dsDNA HS Assay Kit	Thermo Fisher Scientific	Q32851

Before Getting Started

Important Precautions

Read through the entire protocol before beginning library preparation.

- Take note of stopping points throughout the protocol where samples can be safely frozen (-30°C to -10°C) to plan your workflow.
- Use good laboratory practices to prevent contamination of samples by PCR products.
- Use nuclease-free PCR tubes, microcentrifuge tubes and aerosol-barrier pipette tips.
- Wipe down workstation and pipettes with nuclease and nucleic acid cleaning products (such as RNase AWAY, Thermo Fisher Scientific).
- Verify that the thermal cycler used for library preparation is in good working order and currently calibrated according to manufacturer specifications.
- Reaction cleanup with AMPure XP beads is performed at room temperature (20°C to 25°C) and is used repeatedly throughout the Archer library prep workflow. Take care that AMPure XP beads are equilibrated to room temperature and fully resuspended by vortexing until homogenous in both color and appearance prior to drawing out material for each use.

Working with Lyophilized Reaction Pellets

- Archer reagents are provided as individually lyophilized reaction pellets in 0.2mL PCR tube strips.
- Allow pouches to reach room temperature (20°C to 25°C) before opening in order to prevent moisture condensation on tubes.
- Always centrifuge tubes briefly before opening to pull contents down.
- Detach the required number of reaction tubes using a clean razor blade and return any unused portion to the pouch with desiccant packet, reseal and store at 2°C to 8°C.
 - For MBC Adapters and Second PCR tubes remember to label prior to returning unused portions to storage.
- Dissolve, mix and spin down:
 - **Never touch the lysosphere with the pipet tip.**
 - Add sample/reagents to pellet in tubes while on ice.
 - Allow at least 5 seconds for pellets to dissolve.
 - Pipet up and down 8 times to mix after the lysosphere has dissolved.
 - Briefly centrifuge and return to ice before proceeding.

Input Nucleic Acid

- Archer LiquidPlex Library Preparation can tolerate up to 5mM EDTA, however, the optimal starting template is DNA in *EDTA-free* buffer or Ultra-Pure Water.
- This method has been successfully used to produce libraries from 1–300ng total input mass; however, it is advisable to use the maximum allowable input mass (ng) whenever possible. Higher input quantities enable more sensitive variant detection.

- Commercial purification kits recommended by Archer can be found on our website at: <http://archerdx.com/support>.



Reagents to Prepare Before Starting

- Make fresh 10mM Tris-HCl pH 8.0 by mixing 30 μ L 500mM Tris-HCl, pH 8.0 (SA0020) with 1470 μ L Ultra Pure Water (SA0213).
 - 10mM Tris-HCl pH 8.0 is appropriate for use for up to one week after mixing
- Make fresh 70% ethanol by adding 14mL 100% ethanol to the bottle labeled Ultra Pure Water for Ethanol Dilution (SA0022).
 - 70% ethanol is appropriate for use for up to one week after mixing.
 - Tightly close the bottle cap to minimize evaporation when not in use.
- Make fresh 5mM NaOH
 - If working from 1M NaOH, add 5 μ L of 1M NaOH to 995 μ L of Ultra Pure Water to yield 5mM final NaOH
 - If working from 5N, add 10 μ L of 5N NaOH to 990 μ L of Ultra Pure Water to yield 50mM NaOH. Mix well and briefly spin down. Take 100 μ L of 50mM NaOH and combine with 900 μ L of Ultra Pure Water to yield 5mM NaOH. Mix well and briefly spin down.

Thermal Cycler Protocols

- Pre-program your thermal cycler with the following protocols.
- Use the appropriate protocols for specific Archer Assays.
- Verify programming prior to initiating runs.

Complete End Repair	Step	Temperature (°C)	Time (minutes)
	1	25	30
	2	4	Hold

Ligation Step 1	Step	Temperature (°C)	Time (minutes)
	1	37	15
	2	4	Hold

Ligation Step 2	Step	Temperature (°C)	Time (minutes)
	1	22	5
	2	4	Hold

Ligation Elution	Step	Temperature (°C)	Time (minutes)
	1	75	10
	2	4	Hold

First PCR Reaction	Step	Temperature (°C)	Time	Cycles
	1	95	3 minutes	1
	2	95	30 seconds	See Product Insert
	3	Varies – See Product Insert		
	4	72	3 minutes	1
	5	4	Hold	1

Second PCR Reaction	Step	Temperature (°C)	Time	Cycles
	1	95	3 minutes	1
	2	95	30 seconds	See Product Insert
	3	Varies – See Product Insert		
	4	72	3 minutes	1
	5	4	Hold	1

* Note: Use 18 cycles for total sample input mass volumes <9ng. If you regularly experience library yields higher than 200nM you can decrease cycle number as needed.

Molecular Barcoding, Sample Indexing, and Multiplexed Sequencing

Molecule-level barcoding (or unique molecule identifier tagging) and sample-level barcoding (also known as index tagging) are both incorporated during Archer MBC ligation. Molecular barcodes are an integral component of the Archer Analysis software suite, (visit <https://archerdx.com/technology-platform/analysis/> for details). Sample barcodes (i.e. index tags) allow pooled libraries to be sequenced simultaneously thereby enabling maximum sequencing throughput and data demultiplexing during downstream bioinformatics analysis.

Sample Multiplexing

- 1) In order to efficiently utilize the throughput of the MiSeq (or other Illumina sequencing platform) as well as prevent low index diversity within your sequencing run, multiple samples should be sequenced simultaneously. Samples can be identified through a combination of two unique nucleotide sequences (see below for more details), which are subsequently read during the sequencing process. The unique nucleotide sequence is often termed an “index”.
- 2) The Archer library preparation reagents for Illumina utilize a combination of two indices to distinguish between samples. The first index is added during Ligation Step 2 and is embedded in the Archer MBC Adapters for Illumina (p5/i5 index). The second index is added in Second PCR and is embedded in MiSeq Index 1 Primers (p7/i7) within the Second PCR reaction pellets.
- 3) In order to maintain appropriate coverage depth, it is recommended that users determine the maximum number of samples that can be run on a MiSeq flow cell (assuming 12 million reads per run using MiSeq reagents v2). In general, larger panels with more targets will require higher sequencing coverage depth and should be run with fewer samples.

Barcode Diversity

- 4) The Illumina MiSeq will work best when index diversity within a run is high. For example, if eight samples are included in a run, and the user chooses to use only one MBC Adapter paired with eight different Index 1 Primers, the run may fail due to low barcode diversity. In this example it is best to use eight different Archer MBC Adapters paired with eight different Index 1 Primers.
- 5) If using more than 48 MBCs, refer to <http://archerdx.com/support> for adapter compatibility.

Visit <http://archerdx.com/support/faqs> for a list of helpful answers to frequently asked questions or contact us directly at tech@archerdx.com

Step 1: Complete End Repair

- 1) Place an appropriate number of Complete End Repair (SA0205) reaction tubes on ice.
- 2) Adjust purified cfDNA samples to a final volume of **50**µL and combine with Complete End Repair lysosphere on ice.
 - a) Dissolve, mix and spin down (See **Working with Lyophilized Reaction Pellets** in the preceding section, **Before Getting Started**).
 - b) Return tubes to ice
- 3) Transfer reactions to a pre-heated thermal cycler and initiate an incubation using the following program and guidelines:
 - a) Heated lid off



End Repair Incubation Conditions:

Step	Temperature (°C)	Time (minutes)
1	25	30
2	4	Hold

- b) Place samples in the thermal cycler and start program.
- c) After the program has reached 4°C, briefly spin down reactions and place on ice.

Cleanup after End Repair

See Important Precautions section above for guidance on working with AMPure XP beads.

- 1) Completely resuspend AMPure Beads by vortexing.
- 2) Add **3X** volume (**150**µL) of AMPure to each reaction.
- 3) Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
- 4) Incubate for **5** minutes at room temperature (20°C to 25°C).
- 5) Briefly spin down tubes.
- 6) Place tubes on the magnet for **4** minutes **or until beads are fully pelleted** against the tube wall.
- 7) Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. (If the pellet becomes dislodged from the magnet and a portion is drawn into the pipet tip, return contents to the tube and repeat magnet incubation step).
- 8) With tubes still on the magnet, add **200**µL of 70% ethanol.
- 9) Incubate for **30** seconds to allow the bead pellet to fully pellet against side of the tubes.
- 10) Remove and discard the supernatant.
- 11) Repeat steps 8-10 for a total of **two washes** in 70% ethanol.
- 12) After the final wash, use a pipette (≤20µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5** minutes at room temperature with open lids. Take care not to over-dry beads as this will significantly decrease overall recovery (yield) of nucleic acid.
- 13) Elute DNA by resuspending beads in **20**µL 10mM Tris-HCl pH 8.0. Vortex the beads to mix.
- 14) Place tubes back on the magnet for **2** minutes.

Step 2: Ligation Step 1

- 1) Place an appropriate number of Ligation Step 1 (SA0196) reaction tubes on ice.
- 2) Transfer **20**µL of purified DNA from **Reaction cleanup after End Repair** Step 14 into Ligation Step 1 tubes. (It is acceptable for a small amount of AMPure beads to be carried over to this step.)
 - a) Dissolve, mix and spin down.
 - b) Return tubes to ice.
- 3) Transfer reactions to a pre-heated thermal cycler and initiate an incubation using the following program and guidelines:



- a) Use a heated lid ($\geq 100^{\circ}\text{C}$)

Ligation Step 1 Incubation Conditions:

Step	Temperature ($^{\circ}\text{C}$)	Time (minutes)
1	37	15
2	4	Hold

- b) Place samples in the thermal cycler and start program.
- c) After the program has reached 4°C , remove tubes from the temperature block, briefly spin down reactions and place on ice.

Cleanup after Ligation Step 1

See Important Precautions section above for guidance on working with AMPure XP beads.

- 1) Completely resuspend AMPure Beads by vortexing.
- 2) Add **3X** volume (**60**µL) of AMPure to each reaction.
- 3) Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
- 4) Incubate for **5** minutes at room temperature (20°C to 25°C).
- 5) Briefly spin down tubes.
- 6) Place tubes on the magnet for **4** minutes **or until beads are fully pelleted** against the tube wall.
- 7) Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. (If the pellet becomes dislodged from the magnet and a portion is drawn into the pipet tip, return contents to the tube and repeat magnet incubation step).
- 8) With tubes still on the magnet, add **200**µL of 70% ethanol.
- 9) Incubate for **30** seconds to allow the bead pellet to fully pellet against the side of the tubes.
- 10) Remove and discard the supernatant.
- 11) Repeat steps 8-10 for a total of **two washes** in 70% ethanol.
- 12) After the final wash, use a pipette (≤ 20 µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5** minutes at room temperature with open lids. Take care not to over-dry beads as this will significantly decrease overall recovery (yield) of nucleic acid.
- 13) Elute DNA by resuspending beads in **42**µL 10mM Tris-HCl pH 8.0.
- 14) Place tubes back on the magnet for **2** minutes.

Step 3: MBC Adapter Incorporation

- 1) Label MBC Adapter tubes with:
 - a) The sample index tag letter (A, B, or C) and number (1-48) from the MBC Adapters pouch label
 - b) Use a permanent laboratory marker taking care to orient lid hinges to the back as illustrated below:



Important This step incorporates the **P5** index tag for sample-level tracking, take care to record this number for future reference and record which MBC adapter is being used for each sample. Unused tubes must be labeled before returning to the pouch.

- 2) Place an appropriate number of MBC Adapter reaction tubes on ice.
- 3) Combine **40**µL of each DNA sample from **Reaction cleanup after Ligation Step 1** step above with a unique Archer MBC Adapter reaction tube. **Avoid pipetting AMPure beads into this reaction.** If minute amounts of AMPure beads were carried over, simply place MBC Adapter tubes on magnet for one minute and transfer all liquid to the next tubes while MBC Adapter tubes remain on the magnet.
 - a) Dissolve, mix and spin down.
- 4) Immediately proceed to **Step 4: Ligation Step 2.**
- 5) (*Optional step*) Use **1**µL of the remaining DNA sample for fragment analysis.

Step 4: Ligation Step 2

- 1) Place an appropriate number of **Ligation Step 2** (SA0197) reaction tubes on ice.
- 2) Transfer the entire volume of each purified DNA sample from **MBC Adapter Incorporation** Step 3 above to Ligation Step 2 tubes.
 - a) Dissolve, mix and spin down.
 - b) Return tubes to ice.
- 3) Transfer reactions to a pre-heated thermal cycler and initiate an incubation using the following program and guidelines.



- a) Heated lid off

Ligation Step 2 Incubation Conditions:

Step	Temperature (°C)	Time (minutes)
1	22	5
2	4	Hold

- b) Place samples in the thermal cycler and start program.
- c) After the program has reached 4°C, briefly spin down reactions and place on ice.



Safe Stopping Point: It is OK to stop and store the reactions at -30°C to -10°C.

Cleanup after Ligation Step 2

Prepare Ligation Cleanup Beads:

- 1) Completely resuspend Ligation Cleanup Beads by vortexing.
- 2) For each library, pipette **50**µL of Ligation Cleanup Beads (SA0210) into a new 0.2mL PCR tube.
- 3) Place tube(s) on the magnet for **1** minute or until the beads are pelleted.
- 4) Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. (If the pellet becomes dislodged from the magnet and a portion is drawn into the pipet tip, return contents to the tube and repeat magnetic pelleting step.)
- 5) Pipette **50**µL of Ligation Cleanup Buffer (SA0209) into each tube to resuspend beads.

Ligation cleanup procedure:

Caution: Vortex PCR tubes with fingers firmly placed on all lids as the detergent may compromise sealing of lids.

- 1) Transfer **50**µL of Ligation Step 2 reaction into the tubes with Ligation Cleanup Beads and Buffer. Your total volume will be **100**µL.
 - a) Mix samples by vortexing.
 - b) Incubate reactions at room temperature for **5** minutes.
 - c) Mix samples by vortexing.
 - d) Incubate reactions at room temperature for **5** minutes.
 - e) Briefly spin down tubes.
- 2) Place tubes on the magnet for **1** minute **or until beads are fully pelleted** against the tube wall.
- 3) Carefully pipette off and discard supernatant (**100**µL) without disturbing the beads.
- 4) Wash beads **two times** with Ligation Cleanup Buffer (SA0209). For each wash:
 - a) Resuspend beads in **200**µL Ligation Cleanup Buffer by vortexing, briefly spin down, and place back on magnet for **1** minute.
 - b) Once slurry has cleared, discard supernatant.
- 5) Wash beads once with Ultra-Pure Water (SA0213):
 - a. Resuspend beads in **200**µL Ultra-Pure Water by vortexing, briefly spin down and place back on magnet
 - b. Once the slurry has cleared, discard the supernatant.
 - c. Take care to ensure that all of the supernatant has been removed from the beads
- 6) Elute DNA from ligation cleanup beads:
 - a) Resuspend ligation cleanup beads in the volume of 5mM NaOH specified in the **Product Insert (A)** corresponding to your target enrichment panel.
- 7) Transfer reactions to a pre-heated thermal cycler and initiate an incubation using the following program and guidelines.
 - a) Use a heated lid ($\geq 100^{\circ}\text{C}$)

Ligation Elution Incubation Conditions:

Step	Temperature (°C)	Time (minutes)
1	75	10
2	4	Hold



- b) After sample has reached 4°C, briefly spin down and transfer to the magnet
- c) Allow beads to remain on the magnet and fully pellet against the tube wall while you prepare **Step 5: First PCR**.

Step 5: First PCR

- 1) Place an appropriate number of First PCR HS (SA0359) reaction tubes on ice.
 - a) Label tubes by sample number
- 2) Pipette the volume of GSP1 mix specified in the **Product Insert (B)** into each First PCR tube.
 - a) Spin down and return tubes to ice
- 3) Transfer the entire volume of eluate from **Step 4: Ligation Step 2, 7(c)** into appropriately labeled First PCR tube(s).
 - a) Mix and spin down.
 - b) Return tubes to ice.
- 4) Transfer reactions to a thermal cycler and initiate the program specified in the **Product Insert (First PCR Reaction)**.
- 5) When the run has completed, briefly spin down reactions and place on ice.

Cleanup after First PCR

See “**Important Precautions**” section above for guidance on working with AMPure XP beads.

- 1) Completely resuspend AMPure Beads by vortexing.
- 2) Add **1.2X** volume (**48**µL) of AMPure to each reaction.
- 3) Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
- 4) Incubate for **5** minutes at room temperature (20°C to 25°C).
- 5) Briefly spin down tubes.
- 6) Place tubes on the magnet for **4** minutes **or until beads are fully pelleted** against the tube wall.
- 7) Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. (If the pellet becomes dislodged from the magnet and a portion is drawn into the pipet tip, return contents to the tube and repeat magnet incubation step).
- 8) Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - a) Add **200**µL 70% ethanol
 - b) Incubate for 30 seconds at room temperature (20°C to 25°C)
 - c) Carefully remove ethanol and discard
- 9) After the final wash, use a pipette (≤20µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5** minutes at room temperature with open lids. Take care not to over-dry beads as this will significantly decrease overall recovery (yield) of nucleic acid.
- 10) Elute DNA by resuspending beads in the volume of 10mM Tris-HCl pH 8.0 specified in the **Product Insert (C)**.
- 11) Place tubes back on the magnet for 2 minutes.
- 12) Carefully transfer the volume of purified eluate specified in the **Product Insert (D)** to a new 0.2mL PCR tube or proceed directly to **Step 6: Second PCR**.



Safe Stopping Point: It is OK to stop and store the reactions at -30°C to -10°C (or leave in thermal cycler on 4°C hold for no more than 24 hours).

Step 6: Second PCR

- 1) Place an appropriate number of Second PCR–HS (SA0360) reaction tubes on ice.
 - a) **Molecular Barcoding, Sample Indexing & Multiplexed Sequencing** in the **Before Getting Started** section).



Important This step incorporates the **P7** index tag for sample-level tracking, take care to record this number for future reference and record which **MBC** adapter is being used for each sample. Unused tubes must be labeled before returning to the pouch.

Index 1 (P7) sequence table:

Sample Number	Illumina Index 1 P7/i7 Sequence
1	TAAGGCGA
2	CGTACTAG
3	AGGCAGAA
4	TCCTGAGC
5	GGACTCCT
6	TAGGCATG
7	CTCTCTAC
8	CAGAGAGG

- 2) On ice, transfer the entire volume of First PCR eluate into each Second PCR tube.



Be sure to record which P7 index sequence is used with which sample.

- 3) Pipette the volume of GSP2 mix specified in the **Product Insert (E)** into each Second PCR tube.
 - a) Mix and spin down.
 - b) Return tubes to ice.
- 4) Transfer reactions to a thermal cycler and initiate the program specified in the **Product Insert (Second PCR Reaction)**.
- 5) After the program has reached 4°C, briefly spin down reactions and place on ice or immediately proceed to **Cleanup after Second PCR**. It is also acceptable to leave tubes in the thermal cycler overnight at 4°C.

Cleanup after Second PCR

See “Important Precautions” section above for guidance on working with AMPure XP beads.

- 1) Completely resuspend AMPure Beads by vortexing.
- 2) Add **1.2X** volume (**48**µL) of AMPure to each reaction.
- 3) Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
- 4) Incubate for **5** minutes at room temperature (20°C to 25°C).
- 5) Briefly spin down tubes.
- 6) Place tubes on the magnet for **4** minutes **or until beads are fully pelleted** against the tube wall.
- 7) Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. (If the pellet becomes dislodged from the magnet and a portion is drawn into the pipet tip, return contents to the tube and repeat magnet incubation step).
- 8) Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - a) Add **200**µL 70% ethanol
 - b) Incubate for 30 seconds at room temperature (20°C to 25°C)
 - c) Carefully remove ethanol and discard
- 9) After the final wash, use a pipette (≤20µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5** minutes at room temperature with open lids. Take care not to over-dry beads as this will significantly decrease overall recovery (yield) of nucleic acid.
- 10) Elute DNA by resuspending beads in **24**µL 10mM Tris-HCl pH 8.0.
- 11) Place tubes back on the magnet for 2 minutes.
- 12) Transfer the **22**µL of purified solution to a new 0.2mL PCR tube. Stop or proceed directly to **Quantify, Normalize and Sequence**. Be sure to avoid transferring beads to the fresh tube.



Safe Stopping Point: It is OK to stop and store the reactions at -30°C to -10°C.

Quantify, Normalize, and Sequence

Quantify

- 1) Quantify the concentration of each library using the KAPA Universal Library Quantification Kit
 - a) The recommended average fragment length of LiquidPlex libraries for the size-adjustment calculation is 150bp.
 - b) Archer libraries are very concentrated. You will need to dilute libraries 1:20,000-1:250,000 for quantification with KAPA qPCR.

Normalize

After quantification, pool libraries at equimolar concentrations and load the sequencer according to manufacturer instructions. For reference sample sheets and additional recommendations, visit our website at <http://archerdx.com/mbc-adapters> and <http://archerdx.com/support/faqs>.

Sequence – Loading recommendations

- 2) Loading recommendations are provided below. **The final loading concentration must be optimized by each user.**
 - a) For MiSeq, use the read level sequence in the table below.

Level	Read Length
(R1) Read 1	151
(R2) Index Read 1	8
(R3) Index Read 2	8
(R4) Read 2	151

- b) In addition, a reference sample sheet is available for download at: <http://www.archerdx.com/mbc-adapters>. Fill out the sample sheet according to the MiSeq protocol.
 - c) Load sequencing libraries with 5% PhiX as a starting point, prepared as follows:
 - i) Dilute and denature PhiX to 10pM according to the Illumina protocol.



Note: The amount of PhiX depends on the complexity of the final library pool. A higher concentration of PhiX is recommended for libraries prepared from low input masses resulting in low complexity libraries.

- ii) Begin with a **4nM** pool of your barcoded libraries:
 - (1) Combine **10**µL of the 4nM library pool with **10**µL 0.2N NaOH in a 1.5mL microcentrifuge tube. Vortex briefly to mix and incubate for **5** minutes at room temperature.
 - (2) Add **10**µL 200mM Tris pH 7.0 and vortex briefly to mix.
 - (3) Add **970**µL ice-cold Illumina Hyb buffer and vortex briefly to mix. This makes 40pM library.
 - (4) Refer to the table below for amounts of pooled library, PhiX and Hyb Buffer to combine in a new 1.5mL microcentrifuge tube and vortex briefly to mix.
 - (5) Load the entire volume (1.0mL) of the final pool into the appropriate well of the MiSeq cartridge.

MiSeq Loading Guidelines

Desired concentration of final pool (pM)	13	14	15	16	17	18
1) Denature pooled libraries:						
Pooled 4nM libraries (µL)	10	10	10	10	10	10
0.2N NaOH (µL)	10	10	10	10	10	10

2) Incubate for 5 minutes

3) Neutralize and Dilute to 40pM						
Pooled libraries + NaOH from Step 1 (µL)	20	20	20	20	20	20
200mM Tris pH 7.0 (µL)	10	10	10	10	10	10
Hyb Buffer (µL)	970	970	970	970	970	970

4) Dilute library to desired loading concentration (pM)	13	14	15	16	17	18
40pM Libraries from step 2 (µL)	325	350	375	400	425	450
10pM Denatured PhiX (µL)	68	74	79	84	89	95
Hyb Buffer (µL)	607	576	546	516	486	455

5) Load Final Pool into Cartridge						
Final pool (µL)	1000	1000	1000	1000	1000	1000

- d) For **NextSeq**, load sequencing libraries with 20% PhiX as a starting point, prepared as follows:
- i) Dilute and denature PhiX to 20pM according to the Illumina protocol.
 - ii) Begin with a 4nM pool of your barcoded libraries:
 - (1) Combine **10µL** 4nM library pool with **10µL** 0.2N NaOH in a 1.5mL microcentrifuge tube. Vortex briefly to mix and incubate for 5 minutes at ambient room temperature.
 - (2) Add **10µL** 200mM Tris pH 7.0 and vortex briefly to mix.
 - (3) Add **970µL** ice-cold HT1 buffer and vortex briefly to mix. This makes 40pM library.
 - (4) Refer to the table below for amounts of pooled library, PhiX and HT1 Buffer to combine in a new 1.5mL microcentrifuge tube and vortex briefly to mix.
 - (5) Spin down and load the entire volume (1.3mL) of this final pool in 20% PhiX into the appropriate well of the NextSeq cartridge.

NextSeq Loading Guidelines

Desired concentration of final pool (pM)	1.4	1.6	1.8
1) Denature pooled libraries:			
Pooled 4nM libraries (µL)	10	10	10
0.2N NaOH (µL)	10	10	10

2) Incubate for 5 minutes

3) Neutralize and Dilute to 40pM			
Pooled libraries + NaOH from Step 1 (µL)	20	20	20
200mM Tris pH 7.0 (µL)	10	10	10
Hyb Buffer (µL)	970	970	970

4) Dilute library to desired loading concentration (pM)	1.4	1.6	1.8
40pM Libraries from step 2 (µL)	46	52	58
20pM Denatured PhiX (µL)	23	26	29
Hyb Buffer (µL)	1231	1222	1213

5) Load Final Pool into Cartridge			
Final pool (µL)	1300	1300	1300

Analysis

Upon completion of the run, the data should be analyzed using Archer Analysis at <https://analysis.archerdx.com/>. Archer Analysis may be run via either a local software installation or Archer Unlimited. Visit our website at <https://archerdx.com/technology-platform/analysis/> for more information.

Demultiplex NextSeq? libraries according to recommendations in FAQs:
<https://support.archerdx.com/>

© 2020 ArcherDX, Inc. All rights reserved. LiquidPlex™, AMP™, PreSeq®, Archer® and FusionPlex® are trademarks of ArcherDX, Inc., NextSeq® and MiSeq® are registered trademarks of Illumina, Inc. Agencourt®, AMPure® and FormaPure® are registered trademarks of Agencourt Biosciences Corporation, a Beckman Coulter company. SYBR®, Life Technologies™, DynaMag™, Thermo Scientific™, Qubit™ and Maxima™ are registered trademarks of Thermo Fisher Scientific, Inc. KAPA Biosystems® is a registered trademark of KAPA Biosystems, Inc. RNase Away™ is a registered trademark of Molecular Bio-Products, Inc. BioRad®, iTaq™, and SsoAdvanced™ are registered trademarks of Bio-Rad Laboratories, Inc. Qiagen® and QuantiTect® are registered trademarks of Qiagen, Inc.

ArcherDX, Inc.

2477 55th Street, Suite 202
Boulder, CO 80301
303-357-9001

<https://archerdx.com/>