

# Prediction of Liabilities with Chronic Drug Treatment in iPSC Derived Cardiomyocytes using an MEA Platform

Christopher Strock, PhD

## *In vitro* Toxicology & Safety Assessment

### Causes of drug failure in development and in the clinic

- Toxicity remains a leading source of attrition in Drug Dx both pre-clinically but also in clinical development
- The industry has responded to this failure by developing early *in vitro* screen to minimize liabilities.
- Cyprotex has been working over the last 7 years developing *in vitro* models to address tissue specific liabilities

Phase	Preclinical	Preclinical	Phase I-III	Phase I-III	Post Approval
Information	Causes of attrition	Causes of attrition	Causes of attrition	Causes of attrition	Withdrawal
Source	ABPI (2008)	Car (2006)	ABPI (2008)	Olson et al. (2000)	Stevens & Baker (2008)
Sample size	155 CD	88 CD	63 CD	82 CD	47 drugs
Cardiovascular	25%	27%	35%	21%	45%
Hepatotoxicity	15%	8%	29%	21%	32%
CNS	12%	14%	2%	21%	2%
Immunotox	7%	7%	10%	11%	2%
GI	5%	3%	2%	5%	2%
Reprotox	9%	13%	5%	1%	2%
Renal	6%	2%	5%	9%	0%
Carcinogenicity	0%	3%	3%	0%	0%

## Introduction



---

- Late stage failures of drugs in clinical trials have significant costs as well as significant safety risk to patients.
- Identification of liabilities early will save money and allow for prioritization of better compounds.
- Identification of ion channel liabilities will only identify specific acute liabilities for cardiac risk
- Telemerized dog studies are expensive and limit the amount of compounds against a target that can be screened. In vitro testing allows for multiple compounds and scaffolds to be tested to identify a safe compound at affordable prices.
- Other than animal studies, there is no effective way to identify chronic dosing effects
- iPSC cardiomyocytes can remain healthy for more than 2 weeks plated on an MEA plate. Due to the fact that it measures electrical activity without addition of any reagents, extended multi-timepoint assays can be run.

## MEA Based Cardiotoxicity Assay

- MEA measure the electrical activity of entire cells and networks of cells
- Allows measurements of overall cardiac function as well as beat rate, sodium amplitude and QT conduction velocity

<b>Instrument</b>	Maestro 48-well MEA system (Axion BioSystems)
<b>Cell Type</b>	Human iPS cell-derived iCell® cardiomyocytes (Cellular Dynamics International) plated and allowed to fully mature and beat synchronously
<b>Assay Details</b>	Five concentrations in duplicate (dependent on customer requirements) Single time point Additional time points and washout (optional)
<b>Data Delivery</b>	Beat rate and number Field potential duration Amplitude Conduction velocity (optional)



## MEA Based Cardiotoxicity

---

### Physiological based toxicity

---


- MEA can be used to predict toxicity through sub-cytotoxic mechanisms which are the main causes of cardiac liability
  - Electrophysiologic effects driven by ion channel and calcium handling mechanisms
    - Identify hERG liabilities, both direct and trafficking
    - Identify other cardiac ion channel effects
    - Identify liabilities caused by calcium handling. e.g. sofosbuvir and amiodarone causing bradycardia
    - Identify arrhythmias and other liabilities from unknown causes
    - Identify responses to cardiac GPCRs
    - All of these are acute 1 hour responses
- Can run long term experiments
  - Assays can be run over days to weeks because no additional dyes or reagents need to be added
  - Identify delayed responses such as delayed arrhythmias or decreases in Na amplitude over time

## Compounds Tested

---








- Pentamidine – Known hERG trafficker
- Vanoxerine – Multichannel blocker failed for Torsades in clinical trial
- BMS-986094 (INX-08189) – Hepatitis C drug that caused death in clinical trial after multiple weeks
- Compound X – Caused loss of Na amplitude and slope after 24 hours with arrhythmias
- Compound Y – Caused delayed arrhythmias after 6 days with delay in repolarization but no hERG effect

## Pentamidine – hERG Trafficking



# Vanoxerine Case Study – Multichannel Blocker

Treatment	Beat Period	Spike Slope	Spike Amplitude	Field Potential Duration (FPD)	Comments
1 $\mu$ M Vanoxerine	103 $\pm$ 0.7%	33 $\pm$ 5.0%	38 $\pm$ 5.0%	120 $\pm$ 4.4%	
316nM Vanoxerine	115 $\pm$ 2.9%	63 $\pm$ 7.8%	65 $\pm$ 9.3%	157 $\pm$ 7.3%	Arrhythmia
100nM Vanoxerine	110 $\pm$ 0.7%	87 $\pm$ 4.1%	88 $\pm$ 3.2%	145 $\pm$ 1.7%	
31.6nM Vanoxerine	106 $\pm$ 0.7%	95 $\pm$ 3.7%	95 $\pm$ 3.8%	114 $\pm$ 5.2%	
10nM Vanoxerine	102 $\pm$ 0.9%	98 $\pm$ 8.9%	97 $\pm$ 7.6%	107 $\pm$ 2.1%	



Multi channel blocker was in a clinical trial for atrial fibrillation. Dose was increased which caused torsades in two patients. Reached 0.85  $\mu$ M.  
Tested concentrations up to 1  $\mu$ M and timepoints up to 4 hours



## BMS-986094 caused Heart failure in patient

---

### After Patient's Death, Study Shows HCV Drug Cardiotoxic in 14 of 34 Treated Patients

Michael O'Riordan  
September 24, 2014

DURHAM, NC – The development of **BMS-986094** (Bristol-Myers Squibb), a nucleotide-polymerase inhibitor for the treatment of chronic hepatitis C (HCV) infection, was terminated back in August 2012, but a new retrospective review analyzing adverse events provides an in-depth look at the cardiotoxicity associated with the investigational agent.

The findings may have profound implications for other drugs in late-stage testing or others approved for use. Investigators say understanding the mechanisms of cardiac dysfunction and its relevance to other agents will require more study.

Development of BMS-986094 was stopped after a 25-year-old male treated with 200-mg dose experienced rapidly progressive heart failure and died. The patient had been hospitalized with shortness of breath and was found to have a left ventricular ejection fraction (LVEF) <10%.

Now, as part of a retrospective review of the phase 2 study that included the index death, researchers report that 14 of 34 patients treated with BMS-986094 had some evidence of cardiac dysfunction.

The review was published online September 24, 2014 in *Hepatology*.

In total, six patients had severe LV dysfunction (LVEF <30%) and eight had moderate LV dysfunction (LVEF 30%-50%) in one or more evaluations in the six months following drug discontinuation.


- In a phase 2 clinical trial, after 39 days of treatment a patient died from heart failure
- 14 of 34 patients treated at the same time also had some evidence of cardiac dysfunction (LVEF of less than 50%) in the 6 months following discontinuation.
- Chronic treatment with this drug caused unexpected liabilities not previously detected in animal studies
- Could we use human iPSC cardiomyocytes in *in vitro* assays to identify this liability?

## Chronic Dosing of BMS-986094 in MEA assay with Human iPSC Cardiomyocytes

Test Article	Time (hr)	Test Conc. (μM)	% of Vehicle			
			Average Beat Period	Average Na+ Slope	Average Na+ Amplitude	Average FPD
INX-08189	1	10	106 ± 5.4%	78 ± 10.8%	80 ± 11.7%	105 ± 6.2%
		2	106 ± 2.7%	112 ± 9.6%	111 ± 7.3%	107 ± 5.5%
		0.4	105 ± 1.8%	91 ± 23.1%	93 ± 26.2%	99 ± 8.5%
		0.08	106 ± 3.4%	122 ± 18.9%	129 ± 16.1%	105 ± 5.3%
		0.016	106 ± 4.6%	115 ± 16.7%	115 ± 16.3%	107 ± 5.9%
		MEC (μM)	NA	10	10	NA
	120	10	ND	ND	ND	ND
		2	100 ± 4.8%	57 ± 26.0%	41 ± 24.3%	102 ± 8.0%
		0.4	85 ± 2.0%	93 ± 35.1%	96 ± 36.0%	82 ± 10.6%
		0.08	88 ± 3.2%	185 ± 47.5%	203 ± 60.5%	86 ± 8.3%
		0.016	96 ± 3.3%	121 ± 50.7%	125 ± 53.1%	92 ± 4.7%
		MEC (μM)	10	0.08	0.08	10
	288	10	ND	ND	ND	ND
		2	ND	ND	ND	ND
		0.4	70 ± 5.1%	33 ± 13.1%	7 ± 4.6%	ND
		0.08	66 ± 3.0%	130 ± 19.5%	138 ± 24.2%	72 ± 6.5%
		0.016	91 ± 6.4%	107 ± 9.9%	107 ± 12.2%	89 ± 6.9%
		MEC (μM)	0.08	0.08	0.08	0.08
	336	10	ND	ND	ND	ND
		2	ND	ND	ND	ND
		0.4	ND	ND	ND	ND
		0.08	62 ± 5.5%	128 ± 11.7%	136 ± 19.1%	68 ± 7.0%
		0.016	93 ± 6.9%	108 ± 2.9%	110 ± 3.1%	90 ± 5.5%
		MEC (μM)	0.08	0.08	0.08	0.08


Compound causes a progressive loss of activity over time. By 14 days, there is no activity left at a concentration of 400nM and the 80nM concentration has a beat rate almost twice the vehicle treated cells suggesting significant activity at this concentration also.

## BMS-986094 MEA Traces



## Mitochondrial Biogenesis

- Some publications have suggested an effect on mitochondrial biogenesis as the mechanism of toxicity



## BMS-986094 Conclusions

---


- Use of the MEA and CDI iPSC derived cardiomyocytes identified a chronic liability that exists for BMS-986094
- The compound caused a loss of beating at 10, 2, and 0.4 $\mu$ M concentrations by 14 days. This was a chronic issue as there was no effect from this compound at 1 hour.
- There was an effect observed even at 80nM. The cells beat rate increased substantially and the Na amplitude increased significantly. This occurred on two separate plates.
- These chronic effects were observed for this compound while other HepC drugs tested at the same time had no accumulating chronic effects
- The calcium flux measurements basically confirmed what was found in the MEA assay with loss at the top two concentrations and very low level calcium flux at the 0.4  $\mu$ M level. The beat rate for the 80nM concentration was recapitulated in this assay.
- There have been reports of effects on mitochondrial biogenesis for INX-08189. Our data show an upregulation of mitochondrial coded protein at 2.5  $\mu$ M concentration and below. Based on timing we observe and the expression of the mitochondrial proteins, it is likely that this is unrelated to the toxicity observed.
- The data here suggest the value of a chronic in vitro assay to determine cardiac safety. The use of multiple assay platforms and multiple time points would result in a significantly improved safety profile for compounds brought forward into animal and human clinical trials. This would also result in a significant savings for companies.

## Compound X – Decrease in Na Amplitude at 24 hours

			% of Vehicle				
Test Article	Time Point	Test Conc. (μM)	Average Beat Period	Average Na+ Slope	Average Na+ Amplitude	Average FPD	Comment
Compound X	1 hour	150	76 ± 3.1%	62 ± 1.8%	65 ± 1.4%	93 ± 1.3%	
		75	85 ± 0.2%	72 ± 10.0%	76 ± 5.9%	97 ± 1.0%	
		30	92 ± 0.2%	82 ± 12.3%	83 ± 9.9%	99 ± 0.3%	
		15	96 ± 0.1%	95 ± 6.2%	95 ± 4.6%	100 ± 1.4%	
		7.5	98 ± 0.1%	85 ± 15.4%	92 ± 5.4%	101 ± 0.5%	
		MEC (μM)	30	75	75	NA	
	24 hours	150	86 ± 0.5%	12 ± 3.4%	22 ± 2.9%	72 ± 1.0%	
		75	94 ± 0.2%	23 ± 0.0%	30 ± 0.9%	83 ± 0.2%	A,C
		30	132 ± 17.1%	17 ± 0.4%	25 ± 0.9%	TND	A,C
		15	144 ± 21.4%	20 ± 2.3%	29 ± 2.7%	TND	A,C
		7.5	157 ± 1.3%	21 ± 1.8%	29 ± 0.9%	92 ± 6.9%	A,C
		MEC (μM)	150	7.5	7.5	15	


Compound X shows limited activity at 1 hour with increased beat rate at top two or three concentrations and decreased Na amplitude at top 2 concentrations. At 24 hours, all concentrations have greatly decreased Na amplitude and slope and at lower concentrations arrhythmias are occurring. There seems to be some multichannel inhibition with unexplained mechanisms driving changes.

# Compound X MEA Traces





## 144 hour dosing of Compound Y – Delayed unexplained arrhythmic activity



- Traces show the time course of the compound Y. The traces show that the repolarization becomes more delayed over time with the largest increase occurring at 6 days.
- It does not appear to be a hERG effect as we see no decrease in the repolarization peak amplitude. You can see multiple depolarization peaks showing arrhythmias.
- At 10 and 30µM, arrhythmias are observed at 6 days as can be observed in the beat plot.



## Conclusions

---

- Use of an MEA cardiomyocyte assay is an effective way of identifying chronic liabilities.
- Pentamidine, a known hERG trafficking inhibitor, was identified using this assay
- Identified vanoxerine liabilities in extended exposure assay
- The chronic MEA assay was effective at identifying compound X which decreased Na amplitude at 24 hours.
  - May be an effect on cell health that is limiting cellular interactions
- The assay identified a compound which caused delayed arrhythmias (5-6 days) of unknown mechanism in telemerized dogs.
- The assay also identified INX-08189, a compound which caused a death in clinical trial, as a significant liability over the course of 14 days.