

The Hamilton logo, featuring the word "HAMILTON" in a bold, blue, sans-serif font, followed by a circular emblem containing a stylized "H" symbol. The background of the entire page is a warm orange color with a blurred image of a laboratory instrument, likely a liquid chromatography system, and a molecular structure model.

Liquid Chromatography-Mass  
Spectrometry (LC-MS)

## **Analysis in Therapeutic Drug Monitoring**

*'The second eBook of the Hamilton series'*

*Updated - Version 2024*

# Foreword

The efficacy and clinical utility of Therapeutic Drug Monitoring (TDM) relies on the specific, sensitive, and reproducible measurement of analytes. Liquid Chromatography coupled with Mass Spectrometry (LC-MS) has emerged as an attractive alternative to traditional immunoassays. Currently, hundreds of therapeutic drugs are analyzed daily in the clinical setting, using this technique and several *in vitro* diagnostic (IVD)-certified kits have entered the market.

We believe that LC-MS has the potential to become the gold standard for TDM, and therefore we have collaborated with leaders in the field to ensure our automation solutions are compatible with their kits, workflows and instruments.

In this eBook, we discuss the impact of LC-MS in TDM, the drivers behind the need for automation during sample preparation for LC-MS, and the solutions that Hamilton offers to customers working in this field.

This eBook is part of a dedicated campaign for LC-MS analysis in TDM, where we aim to provide our readers with interesting educational resources and additional insights into the way our customers are using Hamilton solutions to accomplish their tasks.

I want to thank our Senior Market Segment Leader Analytics and Chromsystems for their valuable contributions to this eBook. We hope you find the content beneficial.

Your kind feedback is always highly appreciated.



Yours sincerely,

**Dr. Gabriela Boza-Moran**

*Senior Team Leader Field Marketing*

*Hamilton Robotics Direct Sales*

[mboza-moran@hamilton.ch](mailto:mboza-moran@hamilton.ch)

# Table of Contents

Exploring the Role of LC-MS in Therapeutic Drug Monitoring .....	04
LC-MS: Why Do We Need Automation? .....	08
What Can Hamilton Offer? .....	11
<b>Application Note 1:</b> Automated Sample Preparation of More than 100 Drugs for Therapeutic Drug Monitoring .....	15
<b>Application Note 2:</b> Automation Sample Preparation for Immunosuppressant Analysis .....	18



## MassSTAR

CE-IVD solution for LC-MS sample preparation (Chromsystems' kits)



## Dual MassSTAR

RUO solution for LC-MS sample preparation



## easyPunch™ STARlet

RUO solution for sample punching and processing

See page 11



# Exploring the Role of LC-MS in Therapeutic Drug Monitoring

Therapeutic Drug Monitoring (TDM) refers to the repeated measurement of an analyte (a therapeutic drug) in a patient's blood or plasma, with the aim of finding the dosage regimen needed to maintain therapeutic concentrations. This methodology is fundamental in clinical pharmacology, an area that studies the effects of the body on the drugs (pharmacokinetics) and the effects of drugs on the body (pharmacodynamics).<sup>1</sup> TDM is generally required for drugs with a [narrow therapeutic index](#) and/or severe side effects such as immunosuppressants.<sup>2</sup>



For drugs with a narrow therapeutic index, small variations in dose or blood concentration may cause therapeutic failures or adverse drug reactions.

Historically, immunoassays have been the gold standard for the analysis of molecules in clinical chemistry. These tests are based on antibody-antigen binding; however they are subject to several limitations<sup>3</sup> such as: (1) limited specificity, due to cross-reactions; (2) limited reproducibility, due to antibody lot-to-lot variations from the same manufacturer and/or discordance between the results obtained from assays from different manufacturers; and (3) limited sensitivity, due to the low linear dynamic range. Furthermore, since every analyte requires a specific antibody, many do not have a commercially available immunoassay.

Liquid Chromatography coupled with Mass Spectrometry (LC-MS) has emerged as an attractive alternative to traditional immunoassays and other popular enzyme-

subtract tests, as it does not suffer from the same limitations found in immunoassays. LC-MS is highly (1) specific and (2) reproducible, as it identifies compounds based on their unique chemical structure without the need to use antibodies, as well as (3) sensitive, due to its high linear dynamic range. These unique features enable the detection of steroids and other analytes that are difficult to measure at low concentrations. Additionally, LC-MS allows for the simultaneous analysis of several (up to hundreds) compounds from the same sample during the same run.<sup>4,5,6</sup> All these advantages have made LC-MS very attractive in the clinical setting, where thousands of samples (blood, plasma, serum, urine, hair, spittle and milk) are analyzed every year for the assessment of multiple parameters.<sup>7</sup>

The initial adoption of Mass Spectrometry (MS) in the clinical setting was limited by the technical complexities of the MS technique used between the 1980s and 1990s (mainly Gas Chromatography-Mass Spectrometry or GC-MS). The new (soft) ionization techniques, the high performance of the system once coupled with front-end Liquid Chromatography (LC) and a second mass analyzer (MS/MS), gave MS the conditions necessary to be a viable and attractive alternative for routine clinical laboratories.<sup>8</sup>

Initially, the use of LC-MS in TDM was focused on immunosuppressive drugs (also known as immunosuppressants), however its use has now extended to the assessment of anticonvulsants, antibiotics, antifungal drugs, antiviral drugs, antidepressants, anti-cancer drugs and cardiac drugs, among others.<sup>6,11</sup> LC-MS is particularly well-suited for the clinical setting when using tandem mass spectrometry (LC-MS/MS) and performing Selected/

Multiple Reaction Monitoring (SRM/MRM), a technique that selects the analytes to be analyzed and quantified with high precision.<sup>9</sup> LC-MS/MS (particularly triple quadrupole MS) is currently one of the most popular methods for TDM.

Laboratories select clinical tests based on their individual requirements. From a regulatory point of view, it is not uncommon for clinical laboratories to use Laboratory-Developed Tests (LDTs, also known as “in-house tests” or “home brew tests”), which — unlike *in vitro* diagnostics (IVD)-certified tests — do not need US Food and Drug Administration (FDA) clearance or a European CE-IVD mark.

An LDT is an *in vitro* diagnostic test designed, manufactured and used within a single laboratory.<sup>12</sup> Their use in the clinic is expected to be limited in the future. By contrast, an IVD-certified test is an *in vitro* diagnostic test “sold as a complete kit that a laboratory purchases from a manufacturer, and comes with all of the procedures and controls to perform the test”.<sup>13</sup> There are several commercial CE-IVD-certified tests for TDM on the market. However, most automated methods are currently run as LDTs.

## References

1. Hallworth, M., & Watson, I. (2017). Therapeutic Drug Monitoring: Clinical Guide. Fourth Edition. Retrieved 19 June 2020, from [https://www.corelaboratory.abbott/sal/learningGuide/ADD-00061347\\_TDM\\_Learning\\_Guide.pdf](https://www.corelaboratory.abbott/sal/learningGuide/ADD-00061347_TDM_Learning_Guide.pdf)
2. Junaid, T., Wu, X., Thanukrishnan, H., & Venkataramanan, R. (2019). Therapeutic Drug Monitoring. *Clinical Pharmacy Education, Practice and Research*, 425-436. doi: 10.1016/b978-0-12-814276-9.00030-1
3. National Cancer Institute. (2020). NCI Dictionary of Cancer Terms. Retrieved 19 June 2020, from <https://www.cancer.gov/publications/dictionaries/cancer-terms/def/immunoassay>
4. Jannetto, P. (2015). Liquid Chromatography Tandem Mass Spectrometry. Retrieved 19 June 2020, from <https://healthcare-in-europe.com/en/news/liquid-chromatography-tandem-mass-spectrometry.html>
5. Bazydio, L., & Mullins, G. (2016). Technical Comparison of Immunoassay and Mass Spectrometry : June 2016 - MedicalLab Management Magazine. Retrieved 19 June 2020, from <https://www.medlabmag.com/article/1290>
6. van der Gugten, J. (2020). Tandem mass spectrometry in the clinical laboratory: A tutorial overview. *Clinical Mass Spectrometry*, 15, 36-43. doi: 10.1016/j.clinms.2019.09.002
7. Keevil, B., Owen, L., & Adaway, J. (2015). Integrating LC-MS/MS into the Clinical Laboratory. Retrieved 19 June 2020, from [https://www.youtube.com/watch?v=\\_gyV1EGzMx4](https://www.youtube.com/watch?v=_gyV1EGzMx4)
8. Jannetto, P., & Fitzgerald, R. (2016). Effective Use of Mass Spectrometry in the Clinical Laboratory. *Clinical Chemistry*, 62(1), 92-98. doi: 10.1373/clinchem.2015.248146
9. Grebe, S., & Singh, R. (2011). LC-MS/MS in the Clinical Laboratory – Where to From Here?. Retrieved 19 June 2020, from <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3052391/>
10. van der Gugten, J. (2020). Tandem mass spectrometry in the clinical laboratory: A tutorial overview. *Clinical Mass Spectrometry*, 15, 36-43. doi: 10.1016/j.clinms.2019.09.002
11. Garg, U., & Zhang, Y. (2016). Mass Spectrometry in Clinical Laboratory: Applications in Therapeutic Drug Monitoring and Toxicology. *Methods in Molecular Biology*, 1-10. doi: 10.1007/978-1-4939-3252-8\_1
12. FDA. (2018). Laboratory Developed Tests. Retrieved 19 June 2020, from <https://www.fda.gov/medical-devices/vitro-diagnostics/laboratory-developed-tests>
13. Mamuszka, H. (2019). The Neverending LDT vs IVD Debate. Retrieved 19 June 2020, from <https://www.thejournalofprecisionmedicine.com/wp-content/uploads/2019/06/jpm219-Mamuszka.pdf>

# LC-MS FOR THERAPEUTIC DRUG MONITORING IN THE CLINICAL SETTING

Therapeutic Drug Monitoring (TDM) is a clinical practice used to optimize individual dosing regimens by measuring specific drugs at designated intervals to maintain their concentration within a target range.<sup>1,2</sup> The need for accurate, precise, and standardized measurement of drugs presents a major challenge for clinical laboratories and the diagnostics industry.<sup>3</sup>

To be suitable for TDM, a drug should satisfy certain criteria, including:

A narrow therapeutic target range

Significant pharmacokinetic variability

A reasonable relationship between blood concentrations and clinical effects

Established therapeutic range and toxic threshold

Availability of cost-effective drug assays

1


A diagnosis is made and a drug is selected. A dosage schedule is designed to reach a target therapeutic plasma concentration. The drug is administered to the patient.



2


Samples are collected from the patient at repeated intervals and tested in a clinical laboratory using Liquid Chromatography coupled with Mass Spectrometry (LC-MS).

*\*Please note that other analytical methods also exist.*



3

Drug concentration in the patient's blood/plasma is determined. A pharmacokinetic model is applied and clinical judgement is used to adjust the dosage if necessary.




# WHAT IS LIQUID CHROMATOGRAPHY-MASS SPECTROMETRY?

Liquid Chromatography coupled with Mass Spectrometry (LC-MS) is a highly specific, sensitive and reproducible analytical technique that combines the physical separation capabilities of Liquid Chromatography (LC) with the mass analysis capabilities of Mass Spectrometry (MS).<sup>4</sup>


1

Samples are extracted/prepared for analysis and then separated using a High-Performance Liquid Chromatography (HPLC) column.<sup>5</sup>




2

The separated sample species are sprayed into an ion source while transitioning into a gas phase.<sup>5</sup> Samples can be ionized through several methods. The most common for TDM analysis is Electrospray Ionization (ESI)



3

Mass spectrometry separates gas phase ions according to their  $m/z$  (mass to charge ratio) value, using electrical and/or magnetic fields to differentiate ions.<sup>5,6</sup> There are different types of mass analyzers. One of the most common for TDM analysis is the triple quadrupole.



4

As ions cross the mass analyzer, the detector counts them and amplifies their signal.<sup>5</sup>

1. Ghiculesco, R., 2008. Abnormal laboratory results: Therapeutic drug monitoring: which drugs, why, when and how to do it. Australian Prescriber, [online] 31(2), pp.42-44. Available at: <https://www.nps.org.au/australian-prescriber/articles/therapeutic-drug-monitoring-which-drugs-why-when-and-how-to-do-it#article> [Accessed 24 June 2020].


2. Kang, J. and Lee, M., 2009. Overview of Therapeutic Drug Monitoring. The Korean journal of internal medicine, [online] 24(1), p.1. Available at: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2687654/> [Accessed 24 June 2020].

3. Brandhorst, G., Oelerich, M., Maine, G., Taylor, P., Veen, G. and Wallemacq, P., 2012. Liquid Chromatography-Tandem Mass Spectrometry or Automated Immunoassays: What Are the Future Trends in Therapeutic Drug Monitoring?. Clinical Chemistry, [online] 58(5), pp.821-825.

Available at: <https://pdfs.semanticscholar.org/f4e1/d4cdfba4122a7fce8984d6f496394ce2f887.pdf> [Accessed 24 June 2020].

4. Technology Networks. 2020. Liquid Chromatography-Mass Spectrometry (LC-MS) For Oligonucleotide Analysis. [online] Available at: <https://www.technologynetworks.com/analysis/white-papers/liquid-chromatography-mass-spectrometry-lc-ms-for-oligonucleotide-analysis-330452> [Accessed 24 June 2020].

5. CHROMacademy. 2020. Mass Spectrometry Fundamental LC-MS Introduction. [online] Available at: <http://www.ecs.umass.edu/eve/background/methods/chemical/Openlit/Chromacademy%20LCMS%20Intro.pdf> [Accessed 24 June 2020].



# LC-MS: Why Do We Need Automation?

There has been a rapid adoption of LC-MS-based testing in the clinical setting over the past few years. Central laboratories, such as the University Hospital of South Manchester, UK, process between 2,000 and 3,000 samples per day (analyzing several analytes per sample).<sup>1,2</sup> Automation, particularly during the sample preparation step, can significantly increase the throughput of clinical laboratories and provide researchers with more time to focus on more valuable tasks.

Sample preparation is a critical step in LC-MS methodology. Due to the specific technology used for detection, mass spectrometers are very sensitive to the presence of salts and charged ions. Additionally, the detection of low-abundance molecules can be masked by the presence of high abundance ones.<sup>3</sup> The quality and reproducibility of the sample preparation step therefore has a significant influence on the final results. This is particularly critical in the clinical setting, where the results obtained from the tests have a direct influence on the patient's treatment.

There is no standard or universal method for the preparation of samples for LC-MS analysis in the clinical setting. Depending on the type and concentration of the analyte, the matrix, the goal of the assessment and the analytical method used, samples might require a very different treatment (e.g. reduction and alkylation, digestion, desalting, protein-precipitation, solid-phase extraction, etc.).<sup>4,5</sup> In the routine clinical setting, the workflows for each type of measurement need to be standardized to ensure reproducibility. Yet this is often hindered by manual processing errors. Even with the availability of experienced, qualified staff (which is an important source of fixed costs),

inter-operator variability remains. While high-throughput is considered the main selling point of automation, in the clinical practice, the need for high reproducibility is the main reason for transitioning from manual to automated workflows.

Current regulatory changes also require a higher level of compliance for the diagnostic tests performed in the clinical setting (beyond investigational use), including LC-MS-based tests. At present, most clinical laboratories develop and validate their LC-MS assays as LDTs. In the US, laboratories with a Clinical Laboratory Improvement Amendments (CLIA) license are allowed to develop and use their own assays, without the need to obtain an IVD certification.<sup>6</sup> Even though the FDA has expressed the intention to better regulate LDTs, no final guidance has been issued since the initial draft guidance of 2014.<sup>7</sup> It must be highlighted, however, that despite the lack of final guidance, the FDA can, and already does, exercise "enforcement discretion" when there are concerns about the validity of the results provided by an LDT.<sup>6</sup>

In the European Union (EU), and since May 2022, a new regulation (EU 2017/746) came into force which changed the regulatory classification of *In Vitro* Diagnostic (IVD) tests and introduced stricter compliance measures (IVDR).<sup>9</sup> These measures cover the tests, kits, instruments, and users. Automation will help clinical laboratories comply with these new regulations to ensure:

- Traceability
- Standardization of workflows
- Elimination/minimization of human error



- Control of user access
- Workflows for handling errors

Reimbursement policies are expected to be based on the use of IVD tests complying with the new regulations. Many laboratories need to seriously evaluate existing instrument providers on regulatory grounds, as this will impact the complexity of any future LDT validation process.

“Platform selection is one of the most critical choices a diagnostic developer can make because it is effectively a partner for life”.

– Hannah Mamuszka, CEO, Alva 10<sup>6</sup>



## References

1. Keevil, B., Owen, L., & Adaway, J. (2015). Integrating LC-MS/MS into the Clinical Laboratory. Retrieved 19 June 2020, from [https://www.youtube.com/watch?v=\\_gyV1EGzMx4](https://www.youtube.com/watch?v=_gyV1EGzMx4)
2. Garg, U., & Zhang, Y. (2016). Mass Spectrometry in Clinical Laboratory: Applications in Therapeutic Drug Monitoring and Toxicology. *Methods in Molecular Biology*, 1-10. doi: 10.1007/978-1-4939-3252-8\_1
3. Reubsaet, L. (2016). Determination of Very Low-Abundance Diagnostic Proteins in Serum Using Immuno-Capture LC-MS-MS. Retrieved 19 June 2020, from <http://www.chromatographyonline.com/determination-very-low-abundance-diagnostic-proteins-serum-using-immuno-capture-lc-ms-ms?pageID=2>
4. ThermoFisher. (2020). Sample Preparation for Mass Spectrometry | Thermo Fisher Scientific - UK. Retrieved 19 June 2020, from <https://www.thermofisher.com/uk/en/home/life-science/protein-biology/protein-biology-learning-center/protein-biology-resource-library/pierce-protein-methods/sample-preparation-mass-spectrometry.html>
5. van der Gugten, J. (2020). Tandem mass spectrometry in the clinical laboratory: A tutorial overview. *Clinical Mass Spectrometry*, 15, 36-43. doi: 10.1016/j.clinms.2019.09.002
6. Mamuszka, H. (2019). The Neverending LDT vs IVD Debate. Retrieved 19 June 2020, from <https://www.thejournalofprecisionmedicine.com/wp-content/uploads/2019/06/jpm219-Mamuszka.pdf>
7. FDA. (2018). Laboratory Developed Tests. Retrieved 19 June 2020, from <https://www.fda.gov/medical-devices/vitro-diagnostics/laboratory-developed-tests>
8. Institut, J. (2020). Laboratory Developed Tests LDT ~ Game Changer for Medical Labs. Retrieved 19 June 2020, from <https://www.johner-institute.com/articles/regulatory-affairs/and-more/laboratory-developed-tests/>
9. Regulation. (2017). Regulation (Eu) 2017/746 of the European Parliament and of the Council. Retrieved 19 June 2020, from <https://eur-lex.europa.eu/legal-content/EN/TXT/?uri=CELEX:02017R0746-20170505>

# ENSURING COMPLIANCE IN TDM THROUGH AUTOMATION

Liquid Chromatography-Mass Spectrometry (LC-MS) based testing is increasing in popularity for clinical laboratories conducting Therapeutic Drug Monitoring (TDM). In response to an increase in this and other types of clinical tests, European regulations for diagnostic tests performed in the clinic have recently changed (from IVDD to IVDR). Stricter regulation measures are also still in discussion in the US.

Clinical laboratories need to adapt to stricter regulatory control for Laboratory Developed Tests (LDTs), which cover the tests, kits, instruments, and users.

Currently, most clinical laboratories develop and validate their own LC-MS assays.



In 2017, Europe introduced the EU 2017/746, which came with new, harsher regulatory measures. These new regulations are applicable from 2022 onwards.

In the US, laboratories with a Clinical Laboratory Improvement Amendments (CLIA) license can develop and use their assays, without the need to obtain an IVD certification.



In Europe, LDTs can be developed and used in certified labs, when no comparable commercial kit is available. LDTs cannot be distributed commercially.



Automation can be introduced to clinical laboratories to help optimize their testing and increase compliance with these regulations, providing:



Better (end-to-end) traceability



Standardization of workflows



Elimination/minimization of human error



Controlled user access



Robust workflows for handling errors




Increased throughput




Better reproducibility





**Figure 1: MassSTAR Deck Layout (optimized for Chromsystems). The deck layout for the MassTox® TDM Series A kit uses a Reagent Carrier for bottles instead of troughs.**



**Figure 2: STAR CL Deck Layout (optimized for Shimadzu).**




Figure 3: Dual MassSTAR, Deck Layout.




Figure 4: [MPE]<sup>2</sup>® Module. Left: Plate Logistics Module. Right: components of the full system.

microcentrifuges that allow full compatibility with most of the workflows for sample preparation.

For customers looking to process dried blood samples, we offer the [easyPunch™ STARlet \(Figure 6\)](#), a specialized platform for automated sample punching and processing. The easyPunch™ recognizes positively stained samples

(e.g. blood on a white card) and negatively stained samples (e.g., saliva, plasma or urine on indicating cards) and makes punches of 1.2, 2, 3 and 6 mm directly into collection plates. The process is fully traceable.

The final decision about which instrument/ARW to use will depend on the specific needs of the customer.




Figure 5: FlipTube's Solutions.




Figure 6: EasyPunch™ STARlet, Deck Layout.

Independent of the selection, all our solutions ensure:

**Precision**, accurate pipetting achieved through our various proprietary technologies:

1. Compression-Induced O-Ring Expansion (CO-RE®), for tip attachment and positioning
2. Liquid Level Detection (LLD), to detect the exact level of liquids in tubes or plates
3. Anti-Droplet Control (ADC), for correctly pipetting volatile organic solvents

**Reproducibility**, consistency ensured through repeated workflow testing.

**Traceability**, automated barcode verification of samples, reagents, plates and tips, as well as dynamic tracking of each aspiration and dispensation step using our Monitored Air Displacement (MAD) and Total Aspiration and Dispense Monitoring (TADM™) technologies.

## APPLICATION NOTE

# Automated Sample Preparation of More than 100 Drugs for Therapeutic Drug Monitoring

Authors: Maximilian Dietz<sup>1</sup>, Martin Waterstraat<sup>1</sup> and Jan Feldmann<sup>2</sup>

<sup>1</sup> Chromsystems Instruments & Chemicals GmbH, Am Haag 12, 82166 Gräfelfing, Germany

<sup>2</sup> Hamilton Bonaduz AG, Via Crusch 8, 7402 Bonaduz, Switzerland

## Introduction

As each patient's absorption and metabolism of therapeutic drugs may vary, correct dosing is crucial to avoid toxic reactions and to keep the drug at therapeutic levels. Here is where Therapeutic Drug Monitoring (TDM) in serum or plasma patient samples comes into play. The gold standard for this analysis is LC-MS/MS; however, manual sample preparation can be a challenge when dealing with a high number of samples. Here, we describe the use of the Hamilton MassSTAR automation system together with *MassTox*<sup>®</sup> TDM Series A to perform automated sample preparation and analysis by LC-MS/MS of more than 100 drugs.

- Complete CE-IVD-certified workflow
- Easy-to-use



Figure 1. Hamilton MassSTAR

## System Description

The MassSTAR (Figure 1) is based on a Hamilton Microlab<sup>®</sup> STARlet with four channels, a CO-RE<sup>®</sup> Gripper, a Barcode Reader and an integrated centrifuge. The deck consists of carriers for samples, reagents, calibrators and controls, two carriers for Deep-Well Plates (DWP) and pipetting tips, as well as a Hamilton Heater Shaker (HHS) (Figure 2). Per run, up to 96 samples including calibrators and controls can be processed by the system. The method has been optimized to enable best performance.




Figure 2. The Deck Layout of the MassSTAR for *MassTox*<sup>®</sup> TDM Series A

- 1 Carrier for calibrators and controls
- 2 Carriers for sample tubes
- 3 Chromsystems Reagent Carrier
- 4 Carrier for 96-well plates and tips
- 5 Carrier for HHS and counterweight plate
- 6 Centrifuge
- 7 Waste
- 8 Barcode Reader

## Kit Description

The sample preparation is based on *MassTox*<sup>®</sup> TDM Series A from Chromsystems, a modular system that consists of three components: the specific parameter sets, the chromatographic column and the BASIC Kit A. The parameter sets contain individual calibrators, controls and internal standards for the specific drugs. The column itself allows the monitoring of nearly 200 parameters without column switching or change of the mobile phases. In combination with an almost identical sample preparation workflow throughout all parameter sets, this helps to minimize the workload in the laboratory. So far, sample preparation of more than 100 parameters can be automated on the MassSTAR, with more to follow. These include antidepressants, antiepileptic drugs, antimycotics, mycophenolic acid and neuroleptics (Test menu for *MassTox*<sup>®</sup> TDM Series A, see below). For a high data accuracy, nearly all parameters are safeguarded by internal standards. The traceability of results is assured by the use of matrix-based calibration material, which is again traceable to certified reference material.

## Workflow

First, all resources are loaded and the barcode of samples and plates is traced. After homogenization, 50 µl of sample/calibrator/control are transferred to a 96 DWP. 25 µl of Extraction Buffer is added and shaken for 2 min at 600 rpm. 250 µl of prepared Internal Standard Solution are pipetted in each well and shaken for 30 s at 1000 rpm. The precipitate is then separated by centrifugation at 2000 x g for 5 min. The supernatant is transferred to a second DWP and diluted with Dilution Buffer. The dilution factor depends on the parameter set and the specific mass spectrometer.

## Technology

One of the key safety features is the Total Aspiration and Dispense Monitoring (TADM™) of the instrument. Each aspiration and dispensing step in each pipetting channel is monitored in real-time. TADM™ verifies the sample transfer with a traceable digital audit trail and identifies errors such as blood clots in the sample or incorrect volumes that have been transferred. A Graphical User Interface (GUI) guides the user through the loading process along with an inventory check (e.g. barcode matching, checking if the plates are correctly positioned, verifying appropriate liquid levels for the reagents). This means that all reagents and consumables are available in sufficient quantities before the run is started, ensuring that no further user intervention is required, thereby providing additional walk-away time. Descriptive screens monitor the current state of the sample preparation process and show the steps the system will perform next. The GUI also prevents invalid user inputs and protects method files against unauthorized manipulation. Additionally, a tailor-made reagent carrier for **MassTox**® TDM Series enables the loading of assay reagents in their original containers. This means that the reagents do not have to be transferred, saving time, preventing contaminations and a mixing up of reagents, and decreasing the evaporation of the liquid.

## Results

The recovery has been determined for all drugs and is dependent on the analyte, ranging from 81% to 115%.

The repeatability, intermediate precision and reproducibility have been measured to determine the overall method's precision. Samples of at least two different concentration levels were prepared ten times and measured in a single sequence. The variation coefficients ranged between 0.6 and 13.7%, depending on respective analytes and concentration level. Preparation and measurement of samples with at least three different concentration levels on different days (20 x 5 or 10 x 2 x 2 design) resulted in variation coefficients between 0.9 and 18.9%. Furthermore, the reproducibility was determined for at least three concentration levels in three different laboratories, on five different days with a sample preparation that was performed five times (3 x 5 x 5 design). The resulting variation coefficients were between 3.2 and 25.7%. Comprehensive data is available in the instruction manual for each parameter set.

The Lower Limit of Quantification (LLOQ) and Higher Limit of Quantification (ULOQ) have been determined by spiking serum and plasma samples with defined amounts of standard substances for all analytes. For all drugs, LLOQ and ULOQ cover the clinically relevant concentration range required for therapeutic drug monitoring.

## Method comparison

Real patient samples have been collected and analyzed by using **MassTox**® TDM Series A. All samples have been processed manually and on the MassSTAR and compared.

Results show an excellent correlation and comparability between the sample preparation of the MassSTAR and manual processing (Examples see Figure 3).




Figure 3. Passing-Bablok Regression comparing Manual and Automated Sample Preparation for the Four Analytes: Lacosamide, Reboxetine, Amisulpride and Citalopram. The graphs demonstrate a high correlation of data.



## Others

Requirements	Part Number	Provider
MassSTAR	806170; CS806170	Hamilton Bonaduz AG; Chromsystems

## The Complete Automated Testing Menu on the MassSTAR

### MassTox® TDM Series A

#### MassTox® TDM Series A Antidepressants 1/EXTENDED (order no. 92713/XT)

Citalopram, N-Desmethylcitalopram, Duloxetine, Fluoxetine, Desmethylfluoxetine, Fluvoxamine, Mirtazapine, N-Desmethylmirtazapine, Paroxetine, Sertraline, N-Desmethylsertraline, Venlafaxine, O-Desmethylvenlafaxine

#### MassTox® TDM Series A Antidepressants 2/Psychostimulants/EXTENDED (order no. 92715/XT)

Atomoxetine, Bupropion, Erythro-Dihydrobupropion, Threo-Dihydrobupropion, Hydroxybupropion, Clomethiazole, Dosulepin, N-Desmethyldosulepin, Methylphenidate, Mianserin, Milnacipran, Moclobemide, Opipramol, Reboxetine, Ritalinic acid, Tianeptine, Tranylcypromine, Trazodone, Vilazodone, Vortioxetine

#### MassTox® TDM Series A Antiepileptic Drugs All-in-One Method (order no. 92721/XT)

Brivaracetam, Carbamazepine, Carbamazepine-10,11-epoxide, 10-OH-Carbamazepine, 10,11-Dihydroxycarbamazepine, N-Desmethylmesuximide, Ethosuximide, Felbamate, Gabapentin, Lacosamide, Lamotrigine, Levetiracetam, Oxcarbazepine, Perampanel, Phenobarbital, Phenylethylmalonamide (PEMA), Phenytoin, Pregabalin, Primidone, Retigabine, Rufinamide, Stiripentol, Sultiame, Theophylline, Tiagabine, Topiramate, Valproic acid, Vigabatrin, Zonisamide

#### MassTox® TDM Series A Antimycotic Drugs (order no. 92722/XT)

Anidulafungin, Caspofungin, Fluconazole, 5-Flucytosine, Isavuconazole, Itraconazole, Hydroxyitraconazole, Ketoconazole, Micafungin, Posaconazole, Voriconazole

#### MassTox® TDM Series A Mycophenolic Acid (order no. 92716)

Mycophenolic acid, Mycophenolic acid glucuronide

#### MassTox® TDM Series A Neuroleptics 1/EXTENDED (order no. 92712/XT)

Aripiprazole, Dehydroaripiprazole, Clozapine, N-Desmethylclozapine, Haloperidol, Olanzapine, N-Desmethylolanzapine, Quetiapine, Norquetiapine, Risperidone, 9-OH-Risperidone

#### MassTox® TDM Series A Neuroleptics 2/EXTENDED 2 (order no. 92714/XT2)

Amisulpride, Benperidol, Brexpiprazole, Bromperidol, Cariprazine, Chlorpromazine, Chlorprothixene, Flupentixol, Fluphenazine, Guanfacine, lloperidone, Levomepromazine, Loxapine, Lurasidone, Melperone, Perazine, Perphenazine, Pimozide, Pipamperone, Promethazine, Prothipendyl, Sertindole, Sulforidazine, Sulpiride, Thioridazine, Ziprasidone, Zotepine, Zuclopenthixol

These products are not available in all countries. Please contact your local Chromsystems sales representative for availability in your region.

Throughout this document, protected product names may be used without being specifically marked as such.

All rights reserved. Chromsystems, the Chromsystems logo and MassTox® are trademarks of Chromsystems. All other trademarks are the sole property of their respective owners.

© 2022 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. AN-2209-04 –11/2022

# HAMILTON

To find a representative in your area,  
please visit:

[www.hamiltoncompany.com/contact](http://www.hamiltoncompany.com/contact)  
Email: [infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

**United States**  
+1-775-858-3000  
**United Kingdom, Ireland**  
+44 121 272 92 80  
**Brazil**  
+55 11 95914 5000  
**China**  
+86 21 6164 6567

**Denmark, Norway,  
Sweden, Finland**  
+46 8410 27 373  
**France**  
+33 184 008 420  
**Germany, Austria,  
Switzerland**  
+49 89 248 804 804

**Benelux**  
+31 40 209 178 0  
**Italy**  
+39 039 930 06 06  
**Japan**  
+81 3 6435 6850  
**Spain, Portugal**  
+34 930 186 262



## APPLICATION NOTE

# Automated Sample Preparation for Immunosuppressant Analysis

Authors: Maximilian Dietz<sup>1</sup>, Martin Waterstraat<sup>1</sup> and Jan Feldmann<sup>2</sup>

<sup>1</sup> Chromsystems Instruments & Chemicals GmbH, Am Haag 12, 82166 Gräfelfing, Germany

<sup>2</sup> Hamilton Bonaduz AG, Via Crusch 8, 7402 Bonaduz, Switzerland

Therapeutic Drug Monitoring (TDM) requires solutions that comply with IVD requirements and are easy-to-use. Chromsystems, together with Hamilton Robotics, has developed the MassSTAR (Figure 1), a CE-IVD-certified automated solution that will soon (2023) also be in-line with IVDR.

- Complete CE-IVD-certified workflow
- Easy-to-use



Figure 1. Hamilton MassSTAR

## Introduction

Immunosuppressive drugs are used to prevent organ transplant rejection and to treat autoimmune diseases. As each patient's absorption and metabolism of the drugs varies, correct dosing is crucial to avoid toxic reactions while still keeping the drug at therapeutic levels, ensuring the patient's well-being.

TDM of the regularly prescribed immunosuppressive drugs cyclosporin A, everolimus, sirolimus and tacrolimus is common in clinical laboratories. The gold standard for the analysis of these drugs in whole blood is LC-MS/MS.

Whole blood samples tend to form blood clots, which is why their handling can be challenging. Process monitoring and high pipetting precision are crucial to ensure the quality and integrity of the results and compliance with regulatory requirements. However, current customized automation solutions tend to neglect these aspects.

## System Description

The MassSTAR is based on a Hamilton Microlab<sup>®</sup> STARlet with four channels, a CO-RE<sup>®</sup> Gripper, a Barcode Reader and a centrifuge integrated to the right of the system. The deck consists of carriers for samples, one carrier for reagents and one for calibrators and controls. In addition, there are two carriers for Deep-Well Plates (DWP), collection plates and pipetting tips as well as a carrier hosting a Hamilton Heater Shaker (HHS) and plates used as counter balances during centrifugation (Figure 2).




Figure 2. Deck Layout of the Hamilton MassSTAR

- 1 Carrier for calibrators and controls
- 2 Carrier for reagents (troughs)
- 3 Carriers for sample tubes
- 4 Carrier for 96-well plates and tips
- 5 Carrier for HHS and counterweight plate
- 6 Centrifuge
- 7 Waste
- 8 Barcode Reader

The system has a capacity of up to 176 samples per run, including controls and calibrators. The application is based on Hamilton's STAR IVD Software. The method has been optimized to enable best performance.

A user-friendly Graphical User Interface (GUI) makes the system easy-to-use and generates output files ready-to-use for all common LC-MS/MS systems.

## CE-IVD-Certified Workflow

### Kit Description

The Chromsystems ONEMINUTE *MassTox*<sup>®</sup> reagent kit for the analysis of immunosuppressants (Chromsystems PN 93900/1200/DWP) is used on this automation platform. It delivers robust, precise and reproducible results with a run time of approximately one minute per sample. The method is completely validated for the majority of tandem mass spectrometers on the market. Sample preparation is reduced to a simple and effective protein precipitation and online purification step (trap column) and, thus, reduces matrix effects drastically. Isotopically labeled internal standards compensate for all residual matrix effects. The use of multilevel calibrators (6PLUS1<sup>®</sup>) add to result accuracy.



Figure 3. Chromsystems ONEMINUTE *MassTox*<sup>®</sup> Immunosuppressant Reagent Kit

### Workflow

First, all resources are loaded and the barcode of samples and plates are traced. Samples in Sarstedt 2.7 ml, 1.2 ml Monovette, Greiner 3 ml Vacuette or other tube types can be processed. They will then be mixed in a series of elaborate mixing steps to ensure perfect homogeneity before 50 µl of each sample is transferred to a dedicated DWP (Chromsystems PN 93956). Total Aspiration and Dispense Monitoring (TADM<sup>™</sup>) ensures proper sample pipetting. Samples that cannot be pipetted properly, e.g. samples with blood clots or short samples, will be recognized, excluded from further processing and flagged.

After sample transfer, 25 µl of internal standard and 100 µl of extraction reagent are added to the plate and agitated at 1,200 rpm for 2 min on the HHS. Afterward, 250 µl precipitation reagent is added to each sample, followed by another incubation step on the HHS for 3 min at 1,200 rpm. Finally, the system transports the plate to the centrifuge, where it is centrifuged for 3 min at 2,000 x g. The supernatant is then transferred via tip transfer into a collection plate (Chromsystems PN 93058). The samples are then ready to be analyzed by LC-MS/MS. For 96 samples, the entire process takes approximately 100 min.

### Technology

Pipetting whole blood tends to be challenging as blood quickly settles down, reducing the homogeneity of the sample. In addition, blood tends to clot, which can generate problems during pipetting (e.g. by blocking tips and causing insufficient sample to be transferred), and lead to cross-contamination.

To ensure homogeneity of the sample and reduce clots, we use an elaborate mixing procedure, aspirating and dispensing the sample at different heights. This enables maximum homogeneity in the sample and dissolves most blood clots.

Immediately after mixing, we transfer the samples using TADM<sup>™</sup>. TADM<sup>™</sup> monitors the pressure curve of the sample on the pipetting channel and recognizes when a transfer is out of the defined boundaries. An error will lead to a retry, and in the case of a repeated error, to the exclusion of the faulty sample (Figure 4).




Figure 4. TADM<sup>™</sup> Functionalities

## Results

### Repeatability & Within-Laboratory Precision

Performance data on repeatability and within-laboratory precision was determined by measuring three different samples and double processing on ten different days and two runs per day. The procedure is based on CLSI EP05-A3 and corresponds to a 10 x 2 x 2 test design (Table 1).

Table 1. Repeatability & within-Laboratory Precision

Substance	Sample	Mean [ $\mu\text{g/l}$ ]	Repeatability		Within-laboratory precision	
			Coefficient of variation	95% confidence interval	Coefficient of variation	95% confidence interval
Cyclosporin A	Low	52.2	2.7%	2.1–3.9%	3.9%	3.0–5.7%
	Middle	263	3.1%	2.4–4.5%	3.9%	3.1–5.3%
	High	580	2.7%	2.1–3.9%	4.0%	3.0–5.8%
Everolimus	Low	2.17	5.8%	4.4–8.4%	6.0%	4.9–7.8%
	Middle	6.89	3.4%	2.6–4.9%	4.2%	3.4–5.5%
	High	15.8	3.0%	2.3–4.3%	3.4%	2.7–4.4%
Sirolimus	Low	2.80	3.8%	2.9–5.5%	4.9%	4.0–6.4%
	Middle	8.11	4.0%	3.0–5.7%	4.4%	3.6–5.7%
	High	19.1	3.0%	2.3–4.4%	3.3%	2.7–4.3%
Tacrolimus	Low	2.87	4.3%	3.3–6.3%	5.3%	4.3–7.1%
	Middle	9.40	3.4%	2.6–5.0%	3.8%	3.1–4.9%
	High	20.6	2.6%	2.0–3.7%	3.1%	2.5–4.1%

### Reproducibility

The performance data on reproducibility was determined at three sites on the basis of three different samples by 5-fold processing on five different days. The procedure is based on CLSI EP05-A3 and corresponds to a 3 x 5 x 5 test design (Table 2).

Table 2. Reproducibility

Substance	Sample	Mean [ $\mu\text{g/l}$ ]	Reproducibility	
			Coefficient of variation	95% confidence interval
Cyclosporin A	Low	49.9	5.2%	3.9–7.7%
	Middle	258	4.3%	3.5–5.7%
	High	581	5.0%	3.4–9.7%
Everolimus	Low	2.19	12.8%	9.9–18.3%
	Middle	6.74	6.3%	5.2–7.9%
	High	15.5	6.2%	5.1–7.9%
Sirolimus	Low	2.79	7.1%	5.8–9.4%
	Middle	8.04	5.8%	4.8–7.4%
	High	18.8	5.8%	4.8–7.2%
Tacrolimus	Low	2.84	8.0%	6.4–10.6%
	Middle	9.14	5.8%	4.5–8.1%
	High	20.3	4.5%	3.7–5.8%

## MassSTAR vs. Reference Method

The automated method's results with the MassSTAR were compared with a published reference method [1]. The results of the Passing-Bablok analysis show that the sample preparation on MassSTAR is comparable to the reference method. Over the entire concentration range, the slope was between 0.95 to 1.05 for all four analytes, cyclosporin A, everolimus, sirolimus and tacrolimus (Figure 5).




Figure 5: Passing-Bablok Analysis of the Automated Sample Preparation using MassSTAR vs. the Reference Method for Cyclosporin A, Everolimus, Sirolimus and Tacrolimus.

## Others

Requirements	Part Number	Provider
MassSTAR	806170; CS806170	Hamilton Bonaduz AG; Chromsystems
MassTox® Immunosuppressants in whole blood ONEMINUTE test	93900/1200/DWP	Chromsystems

[1] Taibon J, van Rooij M, Schmid R, Singh N, Albrecht E, Anne Wright J, Geletneký C, Schuster C, Mörlein S, Vogeser M, Seger C, Pongratz S, Kobold U. An isotope dilution LC-MS/MS based candidate reference method for the quantification of cyclosporine A, tacrolimus, sirolimus and everolimus in human whole blood. *Clinical Biochemistry* 2020;82:73–84.

Throughout this document, protected product names may be used without being specifically marked as such. All rights reserved. Chromsystems, the Chromsystems logo and MassTox® are trademarks of Chromsystems. All other trademarks are the sole property of their respective owners. © 2022 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. AN-2209-03 –11/2022

# HAMILTON

To find a representative in your area, please visit:

[www.hamiltoncompany.com/contact](http://www.hamiltoncompany.com/contact)  
Email: [infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

**United States**  
+1-775-858-3000  
**United Kingdom, Ireland**  
+44 121 272 92 80  
**Brazil**  
+55 11 95914 5000  
**China**  
+86 21 6164 6567

**Denmark, Norway, Sweden, Finland**  
+46 8410 27 373  
**France**  
+33 184 008 420  
**Germany, Austria, Switzerland**  
+49 89 248 804 804

**Benelux**  
+31 40 209 178 0  
**Italy**  
+39 039 930 06 06  
**Japan**  
+81 3 6435 6850  
**Spain, Portugal**  
+34 930 186 262



# PLATFORMS


## OUR PORTFOLIO

### OF AUTOMATED LIQUID HANDLERS



**Microlab® VANTAGE line**

The next generation of liquid handling automation



**Microlab® STAR<sup>V</sup> Line**

The new benchmark on liquid handling automation



**Microlab® STAR line**

Reliable liquid handling for over 20 years



**Microlab® NIMBUS line**

A compact, multi-channel automated liquid handler



**Microlab® Prep™**

Our smallest footprint liquid handler

# Tools for Analytical Chemistry Analysis

HAMILTON'S SOLUTIONS FOR EVERY NEED




## FlipTubes

For closed reactions



## [MPE]<sup>2</sup>

For solid phase  
extraction with positive  
pressure



## Centrifuge

For spinning samples  
up to 4,200 rpm



## ODC

For incubations between  
+4°C and +99°C



## Punching Module

For card punching



## ApH Module

For pH measurement



## Sealer

For adhesive plate  
sealing with optically  
clear or pierceable foil




## Piercing Tips

For accessing  
samples in various  
sealed labware



## Specialized Tips

For "in-tip" sample  
purification



## Xpress Dialyzer

For "on-plate"  
buffer exchange

Integrations from partners' technologies:





Disclaimer: throughout this eBook, protected product names may be used without being specifically marked as such.

**HAMILTON** 

©Updated 2024 (first version 2020)  
Hamilton Bonaduz AG. All rights reserved.

Web: [hamiltoncompany.com/automated-liquid-handling](https://www.hamiltoncompany.com/automated-liquid-handling)

**Hamilton Company  
United States HQ**  
Tel: +1 775 858 3000

**Hamilton Bonaduz AG  
Europe HQ**  
Tel: +41 58 610 10 10

**United Kingdom &  
Ireland**  
Tel: +44 121 272 92 80

**Spain, Portugal**  
Tel: +34 930 186 262

**France**  
Tel: +33 184 008 420

**Italy**  
Tel: +39 039 930 0606

**Sweden, Norway,  
Denmark & Iceland**  
Tel: +46 8410 273 73

**Benelux**  
Tel: + 31 40 209 178 0

**Germany, Switzerland,  
Austria**  
Tel: +49 89 248 804 804

**China**  
Tel: +86 21 6164 6567