

Mitochondria Staining Protocol

Fixation

- 1) Pre-warm buffers in a 37°C water bath.
- 2) Incubate in fixation buffer for 10 min at 37°C.
- 3) Incubate cells with reducing buffer for 7 min at RT.
- 4) Wash cells 2x with PBS.**

Blocking

- 5) Incubate cells with the blocking buffer for 30 min at RT.

Primary Antibody

- 6) Dilute TOMM20 antibody 1:50 in blocking buffer.
- 7) Incubate cells with the primary antibody solution for 1 h at RT in the dark.***
- 8) Wash cells 3x with PBS.

Secondary Antibody

- 9) Dilute Alexa Fluor® 647-anti-rabbit 1:250 in blocking buffer.
- 10) Incubate cells with the secondary antibody solution for 1 h at RT in the dark.***
- 11) Wash 3x with PBS.

Post Fixation

- 12) Incubate cells with the post fixation buffer for 10 min at RT in the dark.
- 13) Wash 3x with PBS.

*Always use freshly made PFA and NaBH₄ buffers.

**At this stage, you can store the sample at 4°C or proceed with staining.

***Incubate on a seesaw rocker to achieve even staining across the sample.

Buffers

Fixation buffer*

4% PFA with 0.2% glutaraldehyde in PBS.

Reducing buffer*




0.1% w/v. Mix 10 mg of NaBH₄ into 10 mL of PBS.

Blocking buffer

5% BSA and 0.1% Tx100 in PBS.

Post fixation buffer*

4% PFA in PBS.



Images | The outer mitochondrial membrane labeled with the import receptor subunit TOMM20-Alexa Fluor® 647

Top tips

See this ONI blogpost for additional information on how to optimize your immunostaining protocol

Sample holder

- We recommend using the ibidi μ -Slide 8 Well Glass Bottom (1.5H glass coverslip bottom). Use 400 μ l/well for all buffer solutions and 200 μ l/well for all antibody solutions.
- If you don't use the ibidi slides, ensure you do not mount the coverslip with mounting media as this will interfere with the autofocus function of the microscope.

Sample Storage (μ -Slide 8 Well Glass Bottom)

- 1) Wash sample 3x with PBS.
- 2) Fill each well to the top with PBS.
- 3) Place parafilm over the wells to seal in the sample. Ensure there are no air bubbles.
- 4) Place the lid onto the sample firmly.
- 5) Store at sample at 4°C.

Chemicals

Paraformaldehyde (PFA)

TritonX-100 (Tx100) (Sigma, T8787)

Glutaraldehyde (Sigma, G7651- 10ML)

NaBH₄ (Sigma, 71320-25G)

Materials

#1.5H coverslips 170 μ m thickness, 35-75 mm (X dimension) by 15-25 mm (Y dimension)

Multi-well chambers

- ibidi (μ -Slide 8 Well Glass Bottom, No. 1.5) available from [here](#)
- Nunc™ Lab-Tek™ II Chambered Coverglass (ThermoFisher, 154534PK) available from [here](#)

Antibodies

TOMM20 antibody (rabbit) (Thermo Scientific, PA552843)

Alexa Fluor® 647 anti-rabbit, F(ab')₂ (Fisher (Invitrogen A21246), 10236552)

This protocol was adapted from: *A. Jimenez, K. Friedl, C. Leterrier. About samples, giving examples: Optimized Single Molecule Localization Microscopy. Methods, 174 (2020), p. 100-114, 10.1016/j.ymeth.2019.05.008*