

Connecting the LightForce Lab to Your iTero Scanner

This article includes the instructions on how to pair LightForce Orthodontics to your iTero lab destination.

Connecting Your iTero Scanner

In order to send scans directly to LightForce through iTero, you need to configure LightForce Orthodontics as a "preferred lab" in iTero.

If you are using an iTero scanner you can simplify the scan submission process by having LightForce Orthodontics added to your menu.


To establish certified connectivity between your practice and the LightForce lab:

- Contact iTero Customer Advocacy at **(800) 577-8767** or email the 24 hour support at iterosupport@aligntech.com and request that the LightForce lab be added to your "preferred lab" list.
- To complete the registration you will be asked to provide the LightForce company ID number, **49657**.
 - Additional information you may need- scanner serial #, patient ID, office company ID.

After confirmation, you will then see the LightForce lab listed in the "send-to" field when completing the Rx. see instructions below.

1. iTero Scanner Instructions

When creating a New Scan in your iTero Scanner make sure to scan under the Case Type '**iRecord**' and Select the Lab Destination to be '**LightForce Orthodontics**' (See Photo)




The screenshot shows the 'New Scan' interface in the iTero Scanner. The top bar includes a 'New Scan' title, a blue 'R' icon, and a user profile for 'Dr. Test, Test'. Below the top bar, the form displays patient information: 'Doctor: Dr. Test, Test' and 'License: * 12345'. The main form area contains several input fields: 'First Name: *', 'Last Name: *', and 'Chart #'. A red box highlights the 'Case Type: *' dropdown menu, which is currently set to 'iRecord', and the 'Ship To: Lightforce Orthodontics' field.

2. LightForce Portal Instructions

When creating a New Patient with LightForce, select the "Send via iTero" to notify the LightForce team that a scan has been sent through the direct lab connection from the scanner.

Select Patient Scan Source


Select a method of sharing patient scans



Send via iTero™ Send via 3Shape™ Upload Now

Note:

If the patient was scanned under '*Invisalign + iRecord*' or '*Invisalign*' the technician can change the case type to '**iRecord**' from the iTero scanner by selecting 'Add Rx' on the patient's file and most recent scan.




Patient: Smith-Doe, John Dr. Test, Test

Name: Smith-Doe, John
Chart Number:
Last Scan: 05/05/2016
New Scan

ID	Scan Date	Doctor	Case Type	Status
9358961	05/05/2016	Dr. Test, Test	iRecord	Completed

View Rx Viewer **Add Rx** Invisalign Outcome Simulator



New Scan Dr. Test, Test

Doctor: Dr. Test, Test License: * 12345

First Name: * Case Type: * iRecord
Last Name: * Due Date:
Chart #: Ship To: Lightforce Orthodontics

Patient must be scanned within **21 days to change the case type from the iTero scanner**. If the patient was scanned over 21 Days, call iTero support to change the

case type to iRecord. If the patient was scanned under '**Invisalign + iRecord**' you can follow the below instructions to export the scans manually from myitero.com.

- Manually upload iTero scans: Click the link to see instructions >>
Myitero.com Exporting Instructions
- If you do not see the option to change the case type in your iTero scanner, you will need to **call iTero support at (800) 577-8767**, ask support to manually enable the "Add Rx" option back into your iTero scanner.