




Validation stickers

The validation stickers' provide a tool to check if the UV-C emission is sufficient.

Additional to the verification of the sensor in the D25, a disinfection cycle can be validated by using UV validation stickers. This validation can be performed per shift with a recommendation of performing it at least once a week.

How to use the validation stickers in the D25:

1. Place 5 new 'yellow' UV validation stickers in 5 different spots (picture 1) of the D25.
2. Place the roller, with the still unused stickers, back into it's original packaging to avoid day light exposure.
3. Close the lid. One disinfection cycle will be performed by the D25.
4. After 25 seconds the lid will automatically open. Check if the UV validation stickers have turned red (see picture 2)
5. Take out the UV validation stickers and stick them on the validation form (see example of form on next page).
6. On this form fill in the date, time of the test , name of the employees who performed the test, and sign.
7. Keep this form next to the device for administration (see example of form on next page).




Picture 1. (Yellow)



Starting color,
no UV-C exposure




Picture 2. (Red)



Successful
cycle

Notes:

- Always keep the unused stickers in the original packaging to restrain contact with daylight
- Please refer to the stickers' original IFU for detailed information



UVsmart

- ✉ info@uvsmart.nl
- in uvsmart
- 🐦 @UVSmartDelft
- ☎ +31 (0)85 060 98 00

DOC-1212 REV. 0

2021 UV Smart BV. All mentioned trademarks and registered trademarks are the property of UV Smart BV. All rights reserved. No part of this document may be reproduced in any form without written permission from the copyright holder.



UVsmart


Date: _____

Time of test: _____

Name : _____

Signature: _____

Place stickers here:



UVsmart

Date: _____

Time of test: _____

Name : _____

Signature: _____

Place stickers here:

