



The first and only FDA cleared device to assess mild traumatic brain injuries (bleed and concussion) at point-of-care.




BrainScope is an FDA cleared, A.I.-powered, mild TBI decision support triage tool



3 FDA cleared algorithms provide objective biomarkers of brain injury


FDA cleared, ruggedized handheld device



EEG signal processing with no connectivity required

Hand-held, device that uses objective markers of brain bleed and concussion to rapidly assess patient status at point-of-care

EEG acquisition system connects headset to device



FDA cleared disposable EEG headset

Algorithms identify “noise” to assure high quality EEG data
No physician EEG interpretation necessary

Physician ordered, technician operated assessments provide actionable results instantly available at point-of-care

Issues with current mTBI diagnostic tools



Structural – Brain Bleed

Head CT is the Standard of Care

Issues with Head CT

- **91% of results are negative**
- Exposes patients to unnecessary radiation
- 3 to 6+ hour average wait-time
- Average net reimbursement \$153
- Cannot detect concussion

4.0M

head CTs are performed

3.6M

head CTs are negative

Functional - Concussion

No Standard of Care

Issues with Concussion Tests^{1,2}

- Single test can lead to misclassification
- Multiple tests
- Inconsistent usage of tests
- Subjective, easily cheated

up to **50%**

wrong result¹

Current Options:

- ImPACT/other neurocognitive tests
- Eye tracking/blinking
- Balance
- Reading a list of numbers
- Phone apps with questionnaires


1. Raab et al. (2020) *Arch Clin Neuropsychol*.
2. Resch et al (2016) *BMJ Open Sport & Exercise Medicine*.

Outstanding clinical performance

A future new gold standard for mTBI triage



Clinical sensitivity as high or higher than other trusted technologies



- Sensitivity is 99% to ≥ 1 cc of blood
- Similar to D-dimer - high confidence in ruling out unnecessary imaging
- **Outstanding negative predictive value of 98%**

Performance you can trust

In the FDA validation study, false negative patients were considered “**not clinically important,**” did not require neurosurgery, and did not return to the hospital within 96 hours. All were Glasgow Coma Score (GCS) 15.

Clinical Outcomes Study findings reported
100% sensitivity to head CT²

References


1. Hanley, et al. (2017) *Acad Emerg Med*.
2. Naunheim, et al. (2019) *Amer J Emerg Med*.
3. Segal, et al. (2007) *Annals Family Medicine*.
4. Anoop, et al. (2009) *Hematology*.
5. Yiadom, et al. (2017) *J American Heart Association*.
6. Smits, et al. (2010) *Radiology*.
7. Daubert & Jeremias. (2010) *Vascular Health and Risk Management* (Review).

Structural and functional assessments reported at the same time from the same EEG signal inputs



Structural


Structural Injury Classifier (SIC)



The Structural Injury Classifier (SIC) accurately, rapidly and objectively **assesses the likelihood of a brain bleed** ($\geq 1cc$)

Functional

Brain Function Index (BFI)



BrainScope's Brain Function Index (BFI) **assesses brain function impairment and scales with severity;** computed at the same time as the SIC

BrainScope is an All-Taker Test

Includes alcohol-impaired and drug-impaired patients¹

Between 30-50% of people with TBIs are **injured while intoxicated** and about one-third are **under the influence of other drugs.**

Indications for Use

- Ages 18 – 85
- within 72 hours of injury
- Glasgow Coma Score (GCS) 13-15

Setup and testing are simple, quick and intuitive



Prep Patient



5 min



Run Test





5 - 15 min

Start to Finish Results in
15-20 minutes

➔ Quick and easy for operators to use

Fits into any workflow and can be operated
by clinical or non-clinical staff

Real-time feedback for operator
assures high quality recording



BrainScope can significantly improve ED throughput



Hospital ED mild brain injury triage process

without **BrainScope**[®]


- 1 Patient check in
- 2 Triage
- 3 Waiting Room
- 4 Clinician/Physician Assessment
- 5 Wait for CT
- 6 Transport to CT
- 7 CT Scan
- 8 Transport from CT/radiologist reviews CT
- 9 Return to ED, CT review, Clinician Reassessment
- 10 Physician Reassessment
- 11 Discharge

Total time 6+ Hours

with **BrainScope**[®]

- 1 Patient check in
- 2 Triage + **BrainScope**
- 3 Clinician Assessment
- 4 Discharge

Much shorter +
adds Concussion data



BrainScope can have an immediate impact on cost and care



>\$300/patient

Average 3rd party
BrainScope
reimbursement
using CPT/APC codes

\$350/hr.

Opportunity cost
benefit for each
hour of patient
wait-time
avoided

\$\$

Continuum of
care
referrals to
appropriate
follow-up care

**More
patients
seen faster**

Increased
throughput and
productivity in ED

**Improved patient
satisfaction**

Get answers without
radiation
and in less time

**Fit with current
initiatives**

**Choosing
Wisely[®]**

**APPROPRIATE
USE CRITERIA**

MIPS
(Merit-Based Incentive Payment System)

Thank you.



BrainScope[®]

For more information, please contact:

Kent Butler

Vice President of Sales – North America

Kent.Butler@brainscope.com

Cell: 586.855.1647