INSTRUCTIONS FOR USE
STERILE CASPAR DISTRACTION SCREWS

Description:

Caspar Distraction Screws are made from surgical grade stainless steel 1.4441 and are threaded at the distal end with a self-tapping thread. They are available in three lengths as follows; 12mm, 14mm and 16mm.

Application/Intended Use/Performance:

Caspar Distraction Screws are used in cervical spinal surgery procedures to allow access to the intervertebral disc during a discectomy procedure (disc removal).

The Caspar Distraction Screws do not make contact with the central nervous system.

The screws are to be used in conjunction with a Caspar Vertebral Body Distraction Set in the cervical spine anatomy between level C1 and C5, to enhance the exposure of disc space for subsequent decompression or facilitate graft insertion into the disc space for spinal fusion.

To be used on patients with the following condition(s).

- Cervical spine instability
- Degenerative disc disease
- Spinal trauma/fractures
- Spinal tumours

Precautions indicated:

- Anterior cervical discectomy and fusion (ACDF).
- At least two screws are required for each procedure, the selection of which length screw is used is dependent on the condition of the one.

Intended User:

Qualified clinicians typically Neurosurgery and Orthopaedic surgeons with the required specialist knowledge to perform Cervical Spinal Surgery.

Patient Population:

For use on patients from adolescence to 70 years old, subject to the agreement of a qualified clinician and the patient’s bone condition.

Packing & Sterilisation:
Caspar Distraction Screws are supplied sterile single use only. They are packed individually in a box of 10, sterilised by Gamma Irradiation in accordance with EN 11137.

Precautions:

- Must only be used for its intended purpose.
- This device is supplied sterile, intended for single use only.
- Not to be re-sterilised once open.
- Re-use of single use devices present a risk of infection to the patient.
- Do not use if the packaging is damaged or opened. Dispose of immediately if packaging is found to be breached in any way.
- The material is Latex Free.

Storage:

- To be stored under normal ambient conditions.
- Keep away from direct sunlight and must be kept dry.
- Life of the product is limited to the expiry date label.

Disposal:

Can be disposed of without any special instructions and poses no danger if standard medical waste disposal regulations are observed.

	Name des Herstellers
---	----------------------

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	Los-Nummer
---	------------

	Ablaufdatum
---	-------------

	Art der Sterilisation
---	-----------------------

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
 1639

MERCIAN
Surgical Instruments of Excellence

Mercian Surgical Supply Co. Ltd. 10 Topaz Way, Topaz Business Park, Bromsgrove, Worcestershire B61 0GD. UK

Tel: +44 (0) 844 879 1133 Email: info@merciansurgical.com

 Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta

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	Manufacturers name
	Manufacturing date
	Lot number
	Expiry date
	Method of sterilisation
	Indicates single use
	Refer to instructions
	Keep dry
	Keep out of direct sunlight
	Do not use if packaging is damaged
	Consult instructions for use

CE 1639

MERCIAN

Surgical Instruments of Excellence

Mercian Surgical Supply Co. Ltd. 10 Topaz Way, Topaz Business Park, Bromsgrove, Worcestershire B61 0GD. UK
Tel: +44 (0) 844 879 1133 Email: info@merciansurgical.com

 Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013 Malta

Date of issue: **December 2020**Revision number: **5**

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