

DESCRIPTION OF PRODUCT: A pelvic holder to prevent sliding. For chair use only.

Rx ONLY



INDICATIONS FOR USE:

- Patients assessed to be at risk of injury from a fall.
- Patients needing a positioning device for added safety while in a chair.
- Patients who have a tendency to slide down in a chair.

CONTRAINDICATIONS:

- **DO NOT** use on a patient who is or becomes highly aggressive, combative, agitated, or suicidal.
- **DO NOT** use on patients with: ostomy, colostomy, or G-tubes; hernias, severe Cardio Obstructive Pulmonary Disease (COPD); or with post-surgery tubes, incisions, catheters, or monitoring lines. These could be disrupted by a restraint.

ADVERSE REACTIONS

- Severe emotional, psychological, or physical problems may occur: if the applied device is uncomfortable; or if it severely limits movement. If symptoms of these problems ever appear for any reason, get help from a qualified medical authority and find a less restrictive, product or intervention.

POSEY SOFT BELT


REF 4125C Chair use only; 4½"W x 16½" (11 cm x 42 cm) belt pad w/6 foot (2 m) straps

APPLICATION INSTRUCTIONS:

⚠WARNING Make sure patient wears proper undergarments to protect skin.

⚠CAUTION Before use, check device for damage. Discard if you have any questions about patient safety.


1. Lay the device on the chair with the narrow side to the back of the chair.
2. Bring the ends of the connecting straps on the narrow end, down behind the chair, and secure out of the patient's reach (fig. 1).
3. Position the patient as far back in the seat as possible, with the buttocks against the back of the chair.
4. Bring the wide part of the pelvic holder up between the patient's legs.
5. Bring the ends of the connecting straps on the wide end, down between the seat and wheelchair sides at a 45-degree angle (fig. 2).
6. Criss-cross the straps and use quick-release ties to attach the straps to the opposite side kick spurs, out of the patient's reach (fig. 3).
7. If the chair has an adjustable seat, secure straps to a movable part of the chair frame, out of the patient's reach.
8. Check that the straps are secure and will not change position, loosen, or tighten if the patient pulls on them, or if the chair is tilted or adjusted.



⚠WARNING

Heed these warnings to reduce the risk of serious injury or death:

- Monitor skin conditions in the groin area frequently. If the patient slides down or forward, pelvic straps may damage the skin.
- There is a risk of chest compression or suffocation, if the patient's body weight is suspended off the chair seat. Use extreme caution with chair cushions. If a cushion dislodges, straps may loosen and allow the patient to slide off the seat (fig. 4).
- Monitor per facility policy to ensure that the patient cannot slide down, or fall off the chair seat and become suspended (fig. 4).



- **STOP USE AT ONCE:** if the patient has a tendency to slide forward or down in the device; or is able to self-release.

APPLICATION INSTRUCTIONS: (PLASTIC INCONTINENCE SHIELD):

1. Insert the pelvic straps through the shield starting at the widest part and pull the shield up until it rests on the foam pelvic piece.

ADDITIONAL SAFETY AND LAUNDERING INSTRUCTIONS ON OTHER SIDE

