

# SUBMICRON CT FOR THE PHARMACEUTICAL INDUSTRY



## FEATURES

### SUBMICRON SPATIAL RESOLUTION WITHOUT COMPROMISE

- **Voxel size: 325 nm**
- **Quasi parallel-beam geometry**
- **Changeable optical lens magnification**

### UNRIVALED HIGH CONTRAST & SPEED

- **Rotating anode source (1.2 kW)**
- **Sensitive sCMOS-based detector**
- **Adaptive X-ray spectra (Cr, Cu, Mo, W)**

### UNIQUE IMAGING REGIMES

- **Phase-contrast**
- **Dual-target source**
- **Field-of-view extension**

## INTRODUCTION

X-ray submicron computed tomography (submicron CT) is one of the most powerful methods for 3D visualization and inspection of any type of sample or product. This non-destructive method provides sufficient resolution and contrast to evaluate any microstructural features, with the ability to resolve structures even below one micron. Moreover, this method requires minimal/no sample preparation, eliminating complicated tasks such as embedding, coating or thin slicing required with other high-resolution methods. The Rigaku nano3DX represents state-of-the-art laboratory-based nanoscale X-ray imaging. This device, when used with deep learning methods, is an unmatched tool for pharmaceutical applications from R&D to production and inspection.

## INSTRUMENT

The nano3DX is a true X-ray microscope (XRM) with the ability to measure relatively large samples at very high resolution. This is accomplished by using a high-powered rotating anode X-ray source and a high-resolution sCMOS X-ray camera. The rotating anode provides for fast data acquisition and the ability to switch anode materials easily to optimize the data acquisition.




Figure 1: nano3DX particle analysis possibilities.

## CT DATA AND RESULTS

The microstructure of solid dosage forms of pharmaceutical products is a critical factor that impacts disintegration and dissolution rates. As such, microstructure will also play a key role in bioequivalence and therapeutic equivalence. Being able to image the microstructure of a solid dosage form allows optimization of production and formulation procedures to achieve a robust dissolution response. If an out-of-specification dissolution result is later observed, analysis of the solid dosage form's internal structure and microstructure can yield a wealth of insights not accessible through traditional analytical approaches, and help resolve mission critical investigations.

New product development can be a highly time-consuming and expensive task. Using Rigaku's nano3DX, this process can be accelerated, providing immediate feedback on a product's internal structure when any discrepancies between expected and actual attributes can be identified.

nano3DX PROVIDES INSIGHT INTO KEY FACTORS SUCH AS:

### Product Morphological Analysis and Structural Analysis

- **Coatings – thickness analysis**
- **Porosity – pores/porous network analysis**
- **Crystallinity – crystalline phases analysis**
- **Aggregates detection**
- **Distribution of API analysis**
- **Dissolution process evaluation**

### Optimization of Manufacturing Process

#### – Defects Detection

- **Cracks**
- **Coatings thickness homogeneity**
- **Contamination**




Figure 2: Coating's morphological analysis—in the internal structure of the analyzed particle, four specific layers were resolvable: core layer, A= immediate drug layer, B= modified drug layer, C= cellulose acetate layer (ordered according to their thickness).





Figure 3: Porosity analysis—in the internal structure of the analyzed particle, a high number of pores were detected, which is reflected by the total porosity of 7.43 percent. Color coding was used according to distance from particle center: Pores in the core layer are yellow and pores in the outer layers are blue.





Figure 4: Coating thickness homogeneity analysis—the thickness of the outer coating layer (cellulose acetate layer) of the analyzed particle varied from 50 to 80  $\mu\text{m}$ , reflecting low manufacturing quality.





Figure 5: Crystallinity analysis—in the internal structure of the analyzed powder, the crystalline phases were resolvable due to high contrast and resolution. The volume of those phase varied from 500 to 1500  $\mu\text{m}^3$ .

CT analysis performed by



RIGAKU EUROPE SE  
Hugenottenallee 167, Neu-Isenburg, 63263, Germany  
T +49 6102 77999 51 | F +49 6102 77999 99  
rese@rigaku.com | rigaku.com