WELCOME

RIGAKU WEBINAR SERIES X-RAY COMPUTED TOMOGRAPHY FOR MATERIALS SCIENCE *GEOLOGY APPLICATIONS* IS STARTING NOW.





Presenter: Aya Takase

Senior Scientist Rigaku Americas Corporation



Host: Tom McNulty

Senior Vice President Rigaku Americas Corporation





You can send us questions during the presentation. They will be addressed at the end of the presentation.





A recording of this webinar will be available. You will receive an email with a link to it tomorrow.



X-RAY COMPUTED TOMOGRAPHY FOR MATERIALS SCIENCE Geology Applications



2 -2

0 🔲 5

GEOLOGY



ge·ol·o·gy

a science that deals with the history of the earth and its life especially as recorded in rocks (Merriam-Webster)















You will learn:

Keys to imaging rocks
Analysis techniques for rocks
Geology applications



WHAT TO CONSIDER?



Resolution & FOV

Contrast





Resolution & FOV

Contrast







Small FOV, small voxel

Large FOV, large voxel



(15)





WHY NOT ALWAYS USE HIGH-RESOLUTION?





(18)



To keep the file size under control.





Run test scans and analysis.



WHAT HAPPENS IF RESOLUTION IS TOO LOW?













9.4 µm







25.2 µm

9.4 µm







$$25.2\ \mu m$$

9.4 µm







9.1 µm









91 µm



Resolution & FOV

Contrast





Carbonate



Sand stone





Mineral	Density [g/cm ³]
Quartz	2.65
Calcite	2.71
Corundum	4.02
Hematite	7.84

Density contrast is not a problem.



THEN WHAT CAN BE A PROBLEM?



Beam hardening

Partial volume effect







Beam hardening

Partial volume effect







Monochromatic radiation



Polychromatic radiation









Cupping/shading artifact








Softer X-rays \rightarrow High absorption rate

Harder X-rays \rightarrow Low absorption rate





High absorption

Low absorption





WHAT SHOULD WE DO?



Numerical correction

Image correction

Physical correction





WHAT IS OUR GOAL?



Image segmentation

Quantitative analysis





Sandstone





Thresholding (Otsu)



Histogram thresholding

Histographic segmentation

Machine/deep learning segmentation















Threshold depends on the distance from the center.







Original image

Distance map





































Thresholding Original image only

Histographic Original + distance map





Thresholding

Histographic

Deep learning



Segmented solid volume









Thresholding 72.9 vol% Histographic 34.9 vol% Deep learning 34.7 vol%



WHICH IS THE BEST METHOD?



Histogram thresholding

Histographic segmentation

Binary High-contrast

Machine/deep learning segmentation



Beam hardening

Partial volume effect

















Machine learning (Random Forest)







Histogram threshold

Machine learning













Partial volume effect









Histogram threshold



Machine learning



Beam hardening

Partial volume effect







WHAT IF WE CARE ABOUT THE SHAPE?











Pixel segmentation

Surface mesh







Pixel segmentation

Surface mesh



Resolution & FOV

Contrast

Analysis purpose




WHAT CAN WE DO WITH CT FOR GEOLOGICAL SAMPLES?



GEOLOGY APPLICATIONS

- Pore structures
- Cracks and fractures
- Phase textures
- Grain size
- Formation process
- Weathering process
- Meteorites imaging
- Fossils imaging

Rigaku



WHAT CAN WE DO WITH A FOSSIL?



Ammonite fossil









Isolated shell





0.5 cm

Surface mesh



CAN WE QUANTIFY VOIDS IN ROCKS?













Voids: 22.9 vol%



Idaho gray



Sandstones void %

Crab orchard



Liver rock



22.9 vol%

9.0 vol%

9.7 vol%



WHAT ABOUT DIFFERENT PHASES?













Voids 2.0 vol%



WHAT ABOUT GRAIN SIZE?



Sandstone





Object separation



Sandstone





Grain size distribution





Mean volume: 5.5 x $10^7 \ \mu m^3$



CAN WE SEE CRACKS?



Cracked rock













CAN WE ANALYZE A PORE NETWORK?



Packed sand

















Thickness

35 µm



215 µm

















0 µm

Path length





Path length





















Indirect path distance = Tortuosity Direct path distance



5.4



Tortuosity



1.0

Tortuosity



1

5





9



Connectivity



1



https://www.theobjects.com/dragonfly/

DRAGONFLY

FEATURES DE key benefits adv

DEEP LEARNING SHOWCASE advanced processing dragonfly result

LEARN s develop ski SUPPORT get help GET DRAGONFLY achieve real results

DRAGONFLY

From next generation rendering to advanced segmentation... New features to advance your results to new levels

Dragonfly 2020.1

GET STARTED WITH A FREE 30-DAY TRIAL




\bigcirc functions K-Meas. version Discrete rangesmall clip ruit-drevatoral Histogram telatry transform telatoral arranter terrational * ""thresholding chose segment extract parameters chosed erable Profile acale App_ Orthographic App spirit controls minitvision 45" seporting histogram o-mins Boolean post-processing automatically presets 2 filtering Selection and the second secon whole SVD stern Area N/bigck cutors-stad Acolerate Support 🦥 Transform 🔌 estenteds 🚽 machine-learning aligne fusion multiple ntilow Otsu project interactively volumetric Quickly Feret projections sublishing CSV Dynamic large blique water ion registration All and a second a train EX Betadataset Automatically segmen Automatically Automatically Filters^{ML} orientation Filters^{ML} orie Tretard orÿstitch include files international international













Home News Documentation Examples Source Code Publications

Import networks extracted using Maximal Ball, iMorph, 3DMA Rock, SNOW, and more.



Left: Berea sandstone network extracted using the Maximal Ball algorithm.

Right: Toray fibrous electrode extracted using the SNOW algorithm.



Developed by the PMEAL Group at the University of Waterloo







Home News Documentation Examples Source Code Publications

Import networks extracted using Maximal Ball, iMorph, 3DMA Rock, SNOW, and more.



Left: Berea sandstone network extracted using the Maximal Ball algorithm.

Right: Toray fibrous electrode extracted using the SNOW algorithm.



"OpenPNM: A pore network modeling package" Gostick et al. Computing in Science & Engineering. 18(4), p60-74 (2016).



244 µm









9.4 µm

CAN WE SEE FLUID IN PORES/VOIDS?



Packed sand + oil



0 min



15 min

Packed sand + oil





20 min

Packed sand + oil





25 min

Packed sand + oil





Packed sand + oil + brine





GEOLOGY APPLICATIONS

- Pore structures
- Cracks and fractures
- Phase textures
- Grain size
- Formation process
- Weathering process
- Meteorites imaging
- Fossils imaging

Rigaku



USEFUL RESOURCES

- Dragonfly webinars http://orss.ca/ytp2
- OpenPNM

http://openpnm.org/











You just learned: Keys to imaging rocks Analysis techniques for rocks Geology applications



ALL IMAGES WERE COLLECTED ON...



nano3DX CT Lab HX CT Lab GX







To learn more ...



Rigaku.com → Contact





PREVIOUS WEBINARS

www.rigaku.com/en/webinars/ x-ray_ct_introduction









Next on X-ray computed tomography Life Science Applications

October 14th Wednesday 11:00 am PDT / 2:00 pm EDT





https://www.rigaku.com/convention/



Q & A SESSION



Aya Takase

Tom McNulty











We'll follow up with your questions.

Recording will be available tomorrow.

Register for the 7th webinar.



THANK YOU FOR JOINING US SEE YOU NEXT TIME!

