One-handed use Pen-style activation Press-release retraction Safe-trimming groove Disable (Optional)

One-handed use	Pen-style activation	Press-release retraction	Safe-trimming groove	Disable (Optional)
----------------	----------------------	--------------------------	----------------------	--------------------

DEVICE DESCRIPTION: The PenBlade is a scalpel which contains a mechanism to extend and retract a blade between an exposed and a retracted safe position. The PenBlade is fitted with either a 10, 11P, or 15 scalpel blade.

INDICATIONS FOR USE: The PenBlade safety scalpel is indicated for the incision of human skin and tissue, and to cut or trim sutures and catheters up to 6 Fr.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician.

WARNINGS: • The PenBlade is a sterile, single-use only product. Do not re-sterilize and/or reuse. • The scalpel blade is sharp and can cause harm if improperly used. • If the tension button is jammed, do not force it to expose the blade. • Do not cover the PenBlade tip with hand or other part of your body at any time.	PRECAUTIONS: • Store the packaged and sealed device in a dark, dry location. • Do not autoclave or expose the device to organic solvents. • Use the PenBlade before the "Use by" date noted on the package. • Always retract the scalpel blade when you are done using it or setting the PenBlade down.
--	--

Manufactured by: ZIEN Medical Technologies, Inc.
2757 South 300 West, Suite F, Salt Lake City, Utah 84115, USA.
Phone: 1-385-444-2666.

European Representative: Wellkamp Ltd
Suite B, 29 Harley Street, London W1G 9QR, England, United Kingdom
Phone: +44 (20) 32876300; Fax: +44 (20) 76811874

Designed by PenBlade, Inc. Salt Lake City, Utah, USA. Phone: 1-801-903-1113 www.penblade.net

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PB5513 Rev. E - DE

Diagram illustrating the five steps of the Penrite system:

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