

**VTM**

# **VIRAL TRANSPORT MEDIA**

## **ROOM TEMPERATURE STORAGE**

The FDA has issued Emergency Use Authorizations (EUA) for numerous molecular and antigen tests for detection of SARS-CoV-2 (COVID-19). Our Viral Transport Media is listed per FDA guidance, "Enforcement Policy for Viral Transport Media during the Coronavirus 2019 (COVID-19) Public Health Emergency." The VTM is manufactured in the LaunchWorks CDMO facility which is FDA registered, cGMP compliant, and ISO 13485:2016 certified. Call or email us to learn more! 978-338-3045 or [info@LaunchWorksCDMO.com](mailto:info@LaunchWorksCDMO.com)

### **Launchworks VTM Recipe:**

- Fetal Bovine Serum (sterile, heat inactivated, USDA approved)
- Hanks Balanced Salt Solution (no phenol red)
- Gentamicin
- Amphotericin B
- Physiological pH range

### **Media Specifications**

- Microbial Sterility Tested
- DNase and Protease free Tested



# LaunchWorks

Letting You Focus on Science

## Viral Transport Media (VTM)

Product	VTM Filled Tubes	1L VTM Filled Bottles	5L VTM Filled Bottles	Standard Kit	Special Kit
Part Number	LWBVTMCT03	LWBVTMCB1L	LWBVTMCB5L	LWBVTMCK01	LWBVTMCKSP
Description	3 mL of VTM in IATA RATED 12 mL tube	1L Bottles	5L Bottles	Labeled tube in biohazard bag	Any combination of additional parts
Bulk	Case of 50 or 500	One (1) bottle or case of six (6)	One (1) bottle or case of six (6)	Cases of 100 or greater	Cases of 100 or greater
Options	Patient Label	None	None	None	See Below

Storage and Transport Temperature: 2 – 30° Celsius, 12 month storage  
8 – 30° Celsius, 6 month storage

## VTM Specimen Kit

- 1 tube with 3 mLs of sterile LaunchWorks VTM
- Biohazard bag - tube labeled per CDC instructions
- Absorbent pad

### Options:

- Nylon-flocked NP swab
- Barcode or QP-coded labels (up to 3)
- Patient Card
- UN 3373 Return Box

## Swab Specifications

- Naseopharyngeal swab
- Individually packaged
- Nylon flocked
- 155mm, 80 mm breakpoint
- FDA-registered Manufacturer
- Ethylene Oxide Sterlization
- 2 year shelf-life
- CE-marked



# LaunchWorks

A Division of Integreon Global

123 Brimbal Ave.  
Beverly, MA 01915  
Phone: 978-338-3045  
LaunchWorksCDMO.com