



MIMcloud[®] System Requirements

MIMTD-554 • 16 Nov 2021

Minimum

- Windows[®] 10 64-bit
- macOS[®] 10.13/10.14
- Intel[®] Core[™] Duo CPU
- 2GB RAM

Recommended

- Windows 10 64-bit
- macOS 10.15/11.0
- Intel Core i5 CPU
- 4GB+ RAM

Supported Browsers

- Google Chrome[™] (recommended)
- Safari[®]
- Firefox[®]
- Microsoft Edge[™]
- Internet Explorer[™] is not recommended

Firewall/Proxy Access

The following addresses should be whitelisted:

- mim-cloud.appspot.com – uses port 80 and 443
- *.s3.amazonaws.com – uses port 80