

# Nearly invisible Aesthetic Products

Precision engineered and manufactured in the U.S.A.



Epoxy Coated

## ULTRAESTHETIC® ARCHWIRES AND CERAMIC BRACKETS

**Compliment aesthetic brackets with Ultraesthetic™ archwires!**

- Manufactured in both G4™ NiTi and S304 Stainless Steel.
- G&H has developed proprietary manufacturing techniques for coating archwires. The epoxy coating is unique to G&H and the most durable in the industry.



**Available Sizes:**


Round	Rectangular		
.012*	.016x.022	.018x.025	
.014	.016x.027**	.019x.025	
.016	.017x.025	.020x.020	
.018	.018x.018	.021x.025	
.020			

\*Available only in G4™ Nickel Titanium  
\*\*Available only in G4™ Nickel Titanium Europa™II

Actual patient, Heather, wearing G&H Vapor brackets and G4 wires.

### Vapor™ Ceramic Brackets

Vapor brackets are desired by patients for their aesthetics and by doctors for their proven clinical performance, strength and secure bonds.




- Vapor's Unique "Portal" Bonding Base
- Adhesive Flows through Portals for Secure Mechanical Bonds

Aesthetic

## LIGATURE TIES/CHAIN

**Available in 33 colors including clear!**



Clear and Pearl elastomers pictured

Tooth-Colored  
**Stainless Steel Ligature Ties**

Preformed in dead soft temper.

**Long**  
.012 (100/pkg) TCLT012

**Pre-cut**  
.012 (100/pkg) TCSLT012L

**Kobayashi Ties**  
Short  
.014 (100/pkg) TSKH014

### Dyna-Link™ Chain

	Clear	Pearl
<b>Continuous</b>	CDLCL1	CDLPL
<b>Short</b>	SDLCL1	SDLPL
<b>Long</b>	LDLCL1	LDLPL

### Bambino-Ties™

.120 diameter  
**84 Sticks/pkg. (1,008 Ties) BTCL84 BTPL84**

### Versa-Ties™

.120 diameter  
**46 Sticks/pkg. (1,012 Ties) VTCL46 VTPL46**

**33 colors available!**

PEARL	GLOW BLUE	MET. PURPLE	FUCHSIA	YELLOW	GLOW	TURQUOISE	WHITE	BLACK
PEARL BLUE	MET. BLUE	LAVENDER	RUBY RED	MET. GOLD	PEARL GREEN	TRANS. GREEN	MET. SILVER	
LT. SKY BLUE	DEEP SKY BLUE	PINK	MET. RED	ORANGE	NEON GREEN	TOOTH	SMOKE	
ROYAL BLUE	NAVY BLUE	FLUOR. PINK	MAROON	FLUOR. ORANGE	GREEN	CLEAR	GRAY	

# Aesthetic Products



Precision engineered and manufactured in the U.S.A.

**GHOrthodontics.com**

**TOLL-FREE 800-526-1026**

**orders@ghortho.com**

G&H Orthodontics 2165 Earlywood Drive Franklin, IN 46131 U.S.A. 800-526-1026 U.S. & Canada +1 317-346-6655 International