

miniPrevail[®] SL

INTERACTIVE SELF-LIGATING BRACKET SYSTEM

The answer to achieve excellent, predictable treatment results in less time.

An **INTERACTIVE** self ligating bracket which uses a NiTi spring clip that is strong enough to deliver continuous and consistent forces while assuring tooth movement control. When using round archwires (up to .020") allows **PASSIVE** engagement during initial leveling and alignment to reduce friction and improve sliding mechanics. When using rectangular archwires (from .016"x.022"), becomes **ACTIVE** to provide effective torque and rotation control toward the end of treatment.

| Intelligent Design | Patient Comfort |
|---|---|
| <p>COLOR ENHANCED I.D. for easy bracket selection</p> <p>ASSEMBLED USING LASER TECHNOLOGY</p> <p>17-4 SURGICAL-GRADE STAINLESS STEEL for greater strength</p> <p>METAL INJECTION MOLDED (MIM) Constructed of a specialized micro-grain stainless steel MIM powder, ideal for creating small, robust appliances that are inherently stronger and more dimensionally accurate than other similar brackets. Ensures exact tolerances and strength of design</p> <p>COMPOUND CONTOUR BASE for increased bond strength and precise bracket placement</p> <p>MICRO-ETCHED 80 GAUGE MESH ANATOMICAL BONDING BASE matches the anatomical curvature of the tooth for maximum contact and a strong, consistent bond</p> | <p>LOW PROFILE</p> <ul style="list-style-type: none">• reduces occlusal interference• enhances comfort <p>AMPLE TIE WING for elastomeric chain and ties</p> <p>FULL SLOT AREA COVERAGE for total tooth movement control</p> <p>Optional hook on cuspids and bicuspids</p> <p>NITI SPRING CLIP for maximum resistance to permanent deformation, durability, and stable performance - prevents unwanted opening, yet opens easily with an explorer, closes with gentle finger pressure</p> |



BRACKETS | BANDS | TUBES | WIRES | SPRINGS | ELASTOMERICS

miniPrevail®SL Brackets

Available in .018 & .022 Roth* or MBT* prescription

| Roth* Rx Maxillary | Torque | Angle | Rotation | M/D | ID |
|--------------------|--------|-------|----------|-------|----|
| Central | 12° | 5° | 0° | .124" | ● |
| Lateral | 8° | 9° | 0° | .106" | ○ |
| Cuspid | -2° | 11° | 0° | .110" | ● |
| Cuspid Hk | -2° | 11° | 0° | .110" | ● |
| Bicuspid | -7° | 0° | 0° | .110" | ●● |
| Bicuspid Hk | -7° | 0° | 0° | .110" | ● |

| Roth* Rx Mandibular | Torque | Angle | Rotation | M/D | ID |
|---------------------|--------|-------|----------|-------|----|
| Anterior | -1° | 0° | 0° | .098" | ●● |
| Cuspid | -11° | 7° | 0° | .110" | ● |
| Cuspid Hk | -11° | 7° | 0° | .110" | ● |
| 1st Bicuspid | -17° | 0° | 0° | .110" | ●● |
| 1st Bicuspid Hk | -17° | 0° | 0° | .110" | ● |
| 2nd Bicuspid | -22° | 0° | 0° | .110" | ●● |
| 2nd Bicuspid Hk | -22° | 0° | 0° | .110" | ● |

| MBT* Rx Maxillary | Torque | Angle | Rotation | M/D | ID |
|-------------------|--------|-------|----------|-------|----|
| Central | 17° | 4° | 0° | .124" | ● |
| Lateral | 10° | 8° | 0° | .106" | ○ |
| Cuspid | -7° | 8° | 0° | .110" | ● |
| Cuspid Hk | -7° | 8° | 0° | .110" | ● |
| Bicuspid | -7° | 0° | 0° | .110" | ●● |
| Bicuspid Hk | -7° | 0° | 0° | .110" | ● |

| MBT* Rx Mandibular | Torque | Angle | Rotation | M/D | ID |
|--------------------|--------|-------|----------|-------|----|
| Anterior | -6° | 0° | 0° | .098" | ●● |
| Cuspid | -6° | 3° | 0° | .110" | ● |
| Cuspid Hk | -6° | 3° | 0° | .110" | ● |
| 1st Bicuspid | -12° | 2° | 0° | .110" | ● |
| 1st Bicuspid Hk | -12° | 2° | 0° | .110" | ● |
| 2nd Bicuspid | -17° | 2° | 0° | .110" | ● |
| 2nd Bicuspid Hk | -17° | 2° | 0° | .110" | ● |



1 case kits (20 brackets) oriented in a kit box.
Brackets are also available in packages of 5.

G&H prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement by the doctor. MBT is a trademark of 3M Unitek Corporation.

Learn More -
GHOrthodontics.com/miniPrevailSL

Not for sale in USA or Canada

Shop our full line of G&H manufactured products @ GHOrthodontics.com



BRACKETS | BANDS | TUBES | WIRES | SPRINGS | ELASTOMERICS

GHOrthodontics.com TOLL-FREE 800-526-1026 or +1 317-346-6655 orders@ghortho.com