Coil Springs

ORTHODONTICS® Quality Manufacturing. Exceptional Service.



G4™ NiTi CLOSED COIL SPRINGS

Closed Coil Springs are available in 9mm and 12mm passive lengths, each with 5 available force values - Feather Light (50g), Extra Light (100g), Light (150g), Medium (200g), Heavy (250g)10 springs per box.

G&H® offers kits of G4™ NiTi springs designed to allow Orthodontists to choose specific force values based on individual applications.

150 Gram Kit

- Features all necessary CC Springs for 150gm load applications.
- Kit includes 12 boxes of springs, 1 box of measuring gauges, inventory kit box and laminated force chart.



All Kit Components Available Separately

- Individual packages of springs (10 springs per pkg)
- Measuring gauge for precise spring prescription (20 per pkg)

G4™ NiTi OPEN COIL SPRINGS

33	pkg. or convenient 21" spool.	
3	I.D.	FORCE
3	.030	Ex-Light
3	.030	Light
5	.036	Med-Light
5	.030	Medium
3	.036	Medium-Hvy
2	.030	Heavy
5	.036	Ex-Heavy
2	045	Facebow

OPEN COIL G4™

Three 7" lengths per

15mm G4™ OPEN COIL SPRINGS 15mm lengths designed for opening spaces up to 15mm.

10/box FORCE .036 Light Medium .036 036 Heavy

G4™ MOLAR DISTALIZING OPEN COIL **SPRINGS** -STOP WOUND deliver gentle 100 grams of force to

achieve average 1mm to 1.5mm distalization per month without patient cooperation. Three 7" lengths per pack, or 21" spool.

G4™ NiTi TAD CLOSED COIL SPRINGS

The unique universal eyelet design fits over virtually any implant screw and snaps in place locking positively on the screw neck.

Closed Coil TAD Springs are available in 11mm, 13mm, 15mm, and 17mm passive lengths, each with 3 available force values - Light (150g), Medium (200g), and Heavy (250g) 10 springs per hox.

200 Gram Kit

- Weighted kit featuring all necessary TAD Springs for 200gm load applications.
- Kit includes 12 boxes of springs, 1 box of measuring gauges, inventory kit box and laminated force chart.

#KTCCOF200

All Kit Components Available Separately

- Individual packages of springs (10 springs per pkg)
- Measuring gauge for spring determination (20 per pkg)

Nickel Titanium and Stainless Steel Open and Closed

COIL SPRING SPOOLS

S304 STAINLESS STEEL COIL SPRINGS

Available in both open and closed coil configurations. Supplied on 21" spools.





G4™ NITI OPEN AND CLOSED COIL **SPRINGS**

deliver ideal forces for desired tooth movements. .030" I.D. coil diameter supplied on a 21" spool.



TOLL-FREE 800-526-1026

GHOrthodontics.com orders@ghortho.com

Coil Springs





Precision engineered and manufactured by G&H in the U.S.A.