


S&T[®] ACLAND[®] MICRO CLAMPS


'The original and best'



“A” or “V” Pattern?

The “V” pattern is the all-purpose clamp, used on all veins, and almost all arteries. The “A” pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the “A” pattern jaw provides extra security against slippage.

- Light and compact design
- Extremely durable
- A unique gripping surface
- Matt or black finish
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws









Ordering details:

SINGLE CLAMPS - packed in pairs

	Matt code		Black code	Clamp size	Vessel size	Clamping force
	00396 V 00397 A		00462 V 00461 A	B-1	0.4 - 1.0mm	10-14g
	00398 V 00399 A		00464 V 00463 A	B-2	0.6 - 1.5mm	17-23g
	00400 V 00401 A		00466 V 00465 A	B-3	1.0 - 2.25mm	35-44g
	00325		00329	HD-S	1.5 - 3.25mm	60-75g
	00286		00327	RD-S	2.0 - 5.0mm	120-140g











Ordering details:

DOUBLE CLAMPS WITH FRAME - frames are used to tie sutures onto during anastomosis

	Matt code		Black code	Clamp size	Vessel size	Clamping force
	00408 V 00409 A		00474 V 00473 A	ABB-1	0.4 - 1.0mm	10-14g
	00410 V 00411 A		00476 V 00475 A	AAB-2	0.6 - 1.5mm	17-23g
	00412 V 00413 A		00478 V 00477 A	ABB-3	1.0 - 2.25mm	35-44g

Ordering details:

DOUBLE APPROXIMATING CLAMPS - used for end to end anastomosis on small vessels

	Matt code		Black code	Clamp size	Vessel size	Clamping force
	00414 V 00415 A		00480 V 00479 A	B-1	0.4 - 1.0mm	10-14g
	00416 V 00417 A		00482 V 00481 A	B-2	0.6 - 1.5mm	17-23g
	00418 V 00419 A		00484 V 00483 A	B-3	1.0 - 2.25mm	35-44g
	00326		00330	HD-S	1.5 - 3.25mm	60-75g
	00287		00328	RD-S	2.0 - 5.0mm	120-140g

Angled micro vessel clamps are useful for areas where access is restricted and line of sight is hampered.

Ordering details:

SINGLE ANGLED CLAMPS - packed in pairs						
	Matt code		Black code	Clamp size	Vessel size	Clamping force
	00312		00363	DC-1	0.4 - 1.0mm	10-14g
	00245 V 00197 A		00365 V 00313 A	DC-2	0.6 - 1.5mm	17-23g
	00246 V 00198 A		00366 V 00314 A	DC-3	1.0 - 2.25mm	35-44g

Acland Clamp Applicator Forceps

- 14cm/7" length



CAF-4 Forceps for smaller sizes of the Acland Clamps: B-1 to B-3

CAF-5 Forceps for larger sizes of the Acland Clamps: HD-S & RD-S



Acland Clamp Container

MV4400-00



Clamp Box

CB-1

With individual compartments for various clamp sizes