

The Biover® microvascular clamp is a high precision instrument and has been developed under the guidance of international surgeons.

- Highest quality (polycarbonate & stainless steel)
- Single or double clamps, sterile ready for use
- Four sizes available for veins & three sizes for arteries
- Latex free and radiolucent
- Exact pressure every use
- Supplied in boxes of 10
- Each box easily identified by a different colour

VEIN CLAMPS

Single Vein Clamps

Double Approximating Clamps

TKFV-1

Clamping Force
15g/mm²



TKFV-2

Clamping Force
15g/mm²



For vessels between 0.2mm to 1mmØ

TKSV-1

Clamping Force
20g/mm²



TKSV-2

Clamping Force
20g/mm²



For vessels less than 1mmØ

TKMV-1

Clamping Force
30g/mm²



TKMV-2

Clamping Force
30g/mm²



For vessels between 1mm to 2mmØ

TKLV-1

Clamping Force
40g/mm²




TKLV-2

Clamping Force
40g/mm²



For vessels over 2mmØ



ARTERIAL CLAMPS

Single Artery Clamps

Double Approximating Clamps

TKS-1

Clamping Force
40g/mm²



TKS-2


Clamping Force
40g/mm²



For vessels less than 1mmØ

TKM-1

Clamping Force
60g/mm²



TKM-2

Clamping Force
60g/mm²



For vessels between 1mm to 2mmØ

TKL-1

Clamping Force
120g/mm²



TKL-2

Clamping Force
120g/mm²



For vessels between 2mm to 4mmØ

