



SURGEONS OPERATING CHAIRS



Coburg Micro Surgery Stool

Suitable for all micro surgery including: Neuro, Ophthalmic & ENT Surgery.

- Arm rests with swivel ball joints, detachable and autoclavable
- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake
- Comfortable Surgeon's seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61cm to 76cm



GS3016-00S

with Arm Rests + Saddle Seat



GS3016-00R

with Arm Rests + Round Seat

SURGEONS OPERATING CHAIRS



Brake

Foot Pump

Coburg Stool

Suitable for all types of surgery including: Neuro, Ophthalmic & ENT Surgery.

- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake
- Comfortable Surgeon's seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61cm to 76cm



GS3014-00S

Saddle Seat **without** armrest



GS3014-00R

Round Seat **without** armrest

Back rest



Foot plate for height adjustment

The Severn Operating Stool with Back Rest

Suitable for all types of surgery

- Comfortable Surgeon's seat
- Adjustable back rest
- 5 anti-static castors
- Height adjustment with foot plate height range 55cm to 78cm

GS3221-67

Severn Stool with Back Rest

Urology version with Back Rest

- Lower height range 48cm to 67cm

GS3221-57



Foot plate for height adjustment

The Avon Operating Stool **without** Back Rest

Suitable for all types of surgery

- Comfortable Surgeon's seat
- 5 anti-static castors
- Hydraulic height adjustment with foot plate, range 55cm to 78cm

GS3221-65

Avon Stool **without** Back Rest

Urology version **without** Back Rest

- Lower height range 48cm to 67cm

GS3221-55

Mercian Surgical Supply Co. Limited

10 Topaz Business Park, Topaz Way, Bromsgrove, Worcestershire B61 0GD. UK.

www.merciansurgical.com

@MercianSurgical Mercian Surgical

Tel: +44 (0) 844 879 1133 Fax: +44 (0) 844 879 1155 Email: info@merciansurgical.com