

TO ASSIST IN VISUALISING MICRO VESSELS DURING ANASTOMOSIS

Description: Silicon strip with 1mm grid on one side and plain on reverse side, available in blue, green and yellow colours. Packed in boxes of 12 pieces.


Application: For use in micro surgical anastomosis of micro vessels under a microscope to allow the surgeon to better visualise the detail of the vessels, being joined.

The Visibility Material can be used either with the grid side up to give the surgeon a scale of the vessel, alternatively the reverse plain side can also be used.


The material has a none reflective surface eliminating distractions from theatre lighting and has a non-tacky finish to aid mobilisation of the Vessel.



Supplied sterile
in boxes of 12 pieces



1mm grid one side and
matt finish on reverse.
Size: 50mm x 25mm



Ordering details:

COLOUR	ordering code
Yellow	VB-2
Green	VB-3
Blue	VB-4

- X-Ray opaque silicon material
- 1mm grid for accurate measurement of vessels
- Latex free material
- Sterile packed in boxes of 12
- Optimal colours to improve visibility of vessel and reduce eye fatigue