



TIDI Shield®

# Polyethylene Protective Gowns AAMI Level 2

## POLYETHYLENE, OVER-THE-HEAD, PROTECTIVE GOWN, ANSI/AAMI PB70 LEVEL 2 PROTECTION

Made in North America, TIDI's Polyethylene Gowns provide AAMI Level 2 protection for caregivers.

- ANSI/AAMI PB70 Level 2 Protection
- Over-the-head style neck offers ease of donning/doffing
- Open back (back of gown is not protective)
- Thumb-loop wrists prevent sleeve slide
- Waist ties for secure fit
- Universal size

CODE	DESCRIPTION	SIZE	UOM
8576A	Polyethylene, Blue, Over-the-Head, Protective Gown, ANSI/AAMI PB70 Level 2 Protection	Universal	75 EA/CS

## CDC PPE GUIDELINES

CDC recommends using additional infection prevention and control practices during the COVID-19 pandemic, along with standard practices recommended as a part of routine healthcare delivery to all patients. These practices are intended to apply to all patients, not just those with suspected or confirmed SARS-CoV-2 infection\*

\*Centers for Disease Control and Prevention. Interim Infection Prevention and Control Recommendations for Healthcare Personnel During the Coronavirus Disease 2019 (COVID-19) Pandemic. July 15, 2020. <https://www.cdc.gov/coronavirus/2019-ncov/hcp/infection-control-recommendations>.



For more information and to place an order, contact your local TIDI representative or TIDI customer service at [excellence@tidiproductions.com](mailto:excellence@tidiproductions.com) or **800-521-1314**. Learn more at [tidiproductions.com](https://tidiproductions.com).

