



TIDI Shield®

Level 2 Face Masks

TIDI's fluid-resistant Level 2 Face Masks are designed with the end user in mind, offering both protection and comfort. Manufactured in the United States, they're made of soft, breathable material and are designed for long-term wear.



- Features comfortable ear loop design for quick and easy donning and long-term wear
- Made of lightweight, latex-free, fiberglass-free, odorless, and non-irritating material, allowing for easy breathing
- Passes ASTM F2100 Requirements for Performance standards:
 - + BFE at $\geq 98\%$
 - + PFE at $\geq 98\%$
 - + Delta-P < 6.0
 - + ASTM F1862 at 120 mmHg
- Made in the USA



CODE	DESCRIPTION	UNIT QTY.
9030	Fluid-resistant procedure face mask with ear loops, Level 2	50 per box/ 500 per case

DESIGNED TO HELP PROTECT BOTH PATIENTS AND CAREGIVERS

TIDI's line of personal protective equipment—which includes face masks, eyewear, face shields, and gowns—are designed to help protect both caregivers and patients against the spread of airborne diseases and illnesses.

To request samples or place an order, contact your TIDI representative or TIDI Customer Service at **800-521-1314**.

