

PACSonWEB 3.3.400 Release note

Module	Release note
Login	Improved password expiration mail
Image Viewing and Manipulation	Improve auto selection when placing spine labels
Studies	Select multiple instances to move to another study
Multi-Monitor	Client side desktop integration with multi-monitor app
Reports	Correction workflow with audio playback, based on Recognosco speech recognition
DICOM Processor	Don't send mails or SMS for migrated studies
Viewer	FIX: Mouse cursor no longer disappears when doing image manipulations
Image Viewing and Manipulation	Scale ruler added to viewer
MPR	FIX: MPR creation from Local image cache
Local Image Cache	FIX: Local image cache accepts studies with empty accession number
Reports	Automatically add text as template in report
Viewer	Allow to define keyboard shortcuts for external reporting flow
Reports	Integration Recognosco speech recognition (finalization)
Home Reading App	Text blocks can be used in combination with Homereading (phone based speech recognition)
Reports	Support pdf report in the order details tab of the patient history
Viewer	Clear indication for prior studies that are more recent than current study
Download	FIX: Study instance EUID adapted correctly when downloading DICOM
Viewer	FIX: cut-off of the bottom pixels of the viewer in iOS Safari solved
Reports	FIX: text block pop-up now showing on Safari and older Edge versions
Image Viewing and Manipulation	FIX: Window level support for IE11
WCF	FIX: Modality Worklist - Characters before '_' no longer removed from Procedure Code
Management	More intuitive UI for labels in viewer
Reports	Improved robustness of ORU outbound communication
WCF	FIX: C-FIND PATIENT level query needs to be on PATIENT root instead of STUDY root

DICOM Processor	FIX: Viewing and downloading of encapsulated PDF document with empty file name tag
Reports	FIX: Copy of report signature when using copy to clipboard function
Studies	FIX: Resend SMS in Study Administration when source setting "Send SMS on first image" is disabled (Norway)
Worklist	Worklist dashboard moved to worklist tab header when in full screen mode
Image Viewing and Manipulation	Improved usability for Window Level on images with 15 or 16 bit grayscale
Viewer	Loop the RP steps when navigating with shortcut keys * Page-Down on last step --> first step * Page-up on first step --> last step
Image Viewing and Manipulation	Add delta with 180° to angle measurement
Image Viewing and Manipulation	FIX: Use middle mouse button to Window Level on MPR
Multi-Monitor	Add links to open report or order from PACSonWEB in RISONWEB (configuration by Dobco)
WCF	FIX: C-MOVE requests no longer fails if the destination cannot process JPEG Lossless
WCF	FIX: Set access when sending known physician without physician number if received through HL7