# Repeater Modules

# GAS MONITORING & ALARM SYSTEMS

Easy to install wall mounted range of repeater modules which connect to Quantum fixed gas monitors to provide a warning to the user regarding the safety of the monitored area with visual and audible alarms when oxygen concentrations in the air surrounding the monitor fall outside pre-set levels.



### KEY FEATURES

# **CryoPanel 1**

- Large, clear, easy to read, battery powered repeater with LCD display for safe readings of O<sub>2</sub>, CO<sub>2</sub>, CO, H<sub>2</sub> concentrations away from the monitored point.
- > Built in visual and audible alarms.

## AirAlert QBA330

Clear easy to read repeater providing traffic light signals (Red/ Amber/Green) regarding the safety of the monitored area.

## CryoPanel 4, 8 and 12

Centralised monitoring at a single point of gas monitors and other laboratory cryogenic equipment (see separate CryoPanel datasheet for full specification).

## **APPLICATIONS**

Quantum repeater modules allow customisation of each system to suit every application. The range offers a CryoPanel 1 battery powered repeater with LCD display to show  $O_2$ ,  $CO_2$ , CO,  $H_2$  concentrations and provide a warning alarm in a safe location away from danger. The AirAlert QBA330 gives a highly visual Red/Amber/ Green traffic light warning of the safety of the monitored area. CryoPanel 4, 8 and 12 extend gas monitoring at a single panel to include many other types of cryogenic equipment including freezer alarms, fan interfaces, safety valves and lone worker alarms.



# Repeater Modules

# OXYGEN MONITORING & ALARM SYSTEM

## **TECHNICAL SPECIFICATIONS**

| CryoPanel 1          |   |
|----------------------|---|
| Dimensions           | H130 x W105 x D60 mm  |
| Weight               | 450g  |
| Alarm Indications    | ALARM – red flashing LED + audible bleep<br>BATTERY – red constant LED + constant tone<br>DISPLAY – Digital LCD 15–25% oxygen |
| AirAlert QBA330      |   |
| Dimensions           | H130 x W105 x D60 mm  |
| Weight               | 450g  |
| Alarm Indications    | Green (Normal), Amber (Low Alarm), Red (Critical<br>Low Alarm) Audible – 85dBA at 30cm  |
| CryoPanel 4, 8 and 1 | 2   |
| Dimensions           | CryoPanel 4 – H180 x W265 x D82 mm<br>CryoPanel 8 – H523 x W295 x D74 mm<br>CryoPanel 12 – H637 x W295 x D74 mm               |
| Weight               | CryoPanel 4 – 750g<br>CryoPanel 8 – 6kg (Stainless Steel Panel)<br>CryoPanel 12 – 7.2kg (Stainless Steel Panel)               |
| Alarm Indications    | LCD – W130 x H75 mm<br>Audible – 85dBA at 30cm  |
|                      |   |

### SERVICE AND SUPPORT

specifications are available upon request.

| Design & installation         | Quantum can advise on how to integrate your existing equipment and/or design systems to suit specific requirements and maximise safety and reliability.                           |
|-------------------------------|---|
| Maintenance and refurbishment | A maintenance check is advised every 12 months.<br>Quantum offers a choice of maintenance contracts and<br>will refurbish/upgrade any product, regardless of age<br>or condition. |
|                               | s a policy of continuous improvement and we reserve<br>hange specifications without prior notice. Full technical  |



## REPEATER OPTIONS

| PART NO.                 | DESCRIPTION  |
|--------------------------|--|
| CryoPanel 1              | Single channel repeater with LCD display to show oxygen concentration & mimic alarms; includes 10m interconnecting cable.  |
| AirAlert QBA330          | Repeater with trafic light indications and audible alarm.  |
| CryoPanel 4,<br>8 and 12 | Cryogenic central monitoring panel for 1–4 (CryoPanel 4), 1–8 (CryoPanel 8), or 1–12 (CryoPanel 12) QFM330 gas monitors and other laboratory monitoring devices. Standard facility for BMS, safety valve control & extraction systems. |

## **ACCESSORIES**

| PART NO.        | DESCRIPTION  |
|-----------------|--|
| AirAlert QAL200 | Siren/Beacon alarm.  |
| QCA200          | CryoAlert button.  |
| QM24VBAT        | Battery back up for mains powered repeaters.   |
| QAD200/GSM      | Auto dialler (transmits alarm to mobile network).  |
| MSB100          | Mains switching box to provide interface<br>between low voltage relays and mains<br>rated devices. |

**Engineering**Quantum Cryogenics
Unit 13 Langley Business Court

Unit 13 Langley Business Cou Worlds End Newbury Berkshire RG20 8RY ervice

Quantum Cryogenics Unit 1A Ethmar Court Mid Devon Business Par Willand **Head Office** 

Quantum Cryogenics Unit B Cranmere Court Exeter Road Industrial Estate Okehampton Devon EX20 1UF

